

Product Datasheet

BHMT Antibody - BSA Free

NBP1-32443

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-32443

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-32443



NBP1-32443

BHMT Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	0.1M Tris (pH 7), 0.1M Glycine, 10% Glycerol

Product Description


Description	Novus Biologicals Rabbit BHMT Antibody - BSA Free (NBP1-32443) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	635
Gene Symbol	BHMT
Species	Human, Mouse
Immunogen	Recombinant protein encompassing a sequence within the center region of human BHMT. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 10 - 1:500, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000

Images


Western Blot: BHMT Antibody [NBP1-32443] - Sample(30 ug whole cell lysate)A:Hep G2 10% SDS PAGE diluted at 1:1000.



Immunocytochemistry/Immunofluorescence: BHMT Antibody [NBP1-32443] - Analysis of methanol-fixed A431, using antibody at 1:500 dilution.



Immunohistochemistry-Paraffin: BHMT Antibody [NBP1-32443] - Paraffin-embedded Colon Carcinoma; antibody dilution 1:500.




Western Blot: BHMT Antibody [NBP1-32443] - Sample (30 ug of whole cell lysate)

A: whole zebrafish

10% SDS PAGE


NBP1-32443 diluted at 1:1000



Immunohistochemistry: BHMT Antibody [NBP1-32443] - Immunohistochemical analysis of paraffin-embedded zebrafish tissue, using BHMT antibody (NBP1-32443) at 1:300 dilution.



Western Blot: BHMT Antibody [NBP1-32443] - Various tissue extracts (50 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with BHMT antibody (NBP1-32443) diluted at 1:1000.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-32443

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBC1-18523	Recombinant Human BHMT His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-32443

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

