

Product Datasheet

TAF15 Antibody - BSA Free

NBP1-32461

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-32461

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-32461



NBP1-32461

TAF15 Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	0.1M Tris (pH 7), 0.1M Glycine, 10% Glycerol
Target Molecular Weight	62 kDa

Product Description


Description	Novus Biologicals Rabbit TAF15 Antibody - BSA Free (NBP1-32461) is a polyclonal antibody validated for use in IHC, WB, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	8148
Gene Symbol	TAF15
Species	Human, Mouse
Reactivity Notes	Negative Xenopus reactivity reported by customer review.
Immunogen	Recombinant protein encompassing a sequence within the center region of human TAF15. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunoprecipitation 1:10 - 1:500, Immunohistochemistry-Paraffin 1:100-1:1000

Images


Western Blot: TAF15 Antibody [NBP1-32461] - (30 ug whole cell lysate)
 A:Raji 7.5% SDS PAGE diluted at 1:1000




Immunocytochemistry/Immunofluorescence: TAF15 Antibody [NBP1-32461] - Paraformaldehyde-fixed HeLa, using antibody at 1:200 dilution.




Immunohistochemistry-Paraffin: TAF15 Antibody [NBP1-32461] - Paraffin-embedded human breast carcinoma. TAFII68 antibody [N1N3] diluted at 1:250.



Western Blot: TAF15 Antibody [NBP1-32461] - (30 ug of whole cell lysate) A:NIH-3T3 7.5% SDS PAGE diluted at 1:1000



Immunohistochemistry-Paraffin: TAF15 Antibody [NBP1-32461] - Paraffin-embedded CL1-0 xenograft, using antibody at 1:100 dilution.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-32461

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-22848	Recombinant Human TAF15 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-32461

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

