

Product Datasheet

HN1 Antibody - BSA Free NBP1-32510

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-32510

Updated 9/25/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-32510



NBP1-32510

HN1 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS, 20% Glycerol
Target Molecular Weight	16 kDa

Product Description	
Description	Novus Biologicals Knockout (KO) Validated Rabbit HN1 Antibody - BSA Free (NBP1-32510) is a polyclonal antibody validated for use in IHC, WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	51155
Gene Symbol	JPT1
Species	Human
Reactivity Notes	Mouse (88%), Rat (82%), Bovine (88%), Dog (83%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human HN1. The exact sequence is proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunohistochemistry, Knockout Validated
Recommended Dilutions	Western Blot 1:500-1:3000, ELISA 1:1000-1:10000, Immunohistochemistry, Immunohistochemistry-Paraffin, Knockout Validated
Application Notes	IHC-P Assay dependent

Images


Western Blot: HN1 Antibody [NBP1-32510] - Wild-type (WT) and HN1 knockout (KO) 293T cell extracts (30 ug) were separated by 12% SDS-PAGE, and the membrane was blotted with HN1 antibody diluted at 1:1000. HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




ELISA: HN1 Antibody [NBP1-32510] - An ELISA plate is coated with 50 uL of HN1 recombinant protein at concentration ranged from 0.004 ug/mL to 1 ug/mL. The coated protein is detected with HN1 antibody at 0.25 ug/mL.




ELISA: HN1 Antibody [NBP1-32510] - An ELISA plate is coated with 50 uL of HN1 recombinant protein at concentration ranged from 0.004 ug/mL to 1 ug/mL. The coated protein is detected with HN1 antibody at 0.25 ug/mL.




Western Blot: HN1 Antibody [NBP1-32510] - Various whole cell extracts (30 ug) were separated by 15% SDS-PAGE, and the membrane was blotted with HN1 antibody (NBP1-32510) diluted at 1:3000.



Western Blot: Rabbit Polyclonal HN1 Antibody [NBP1-32510] - Various whole cell extracts (30 ug) were separated by 15% SDS-PAGE, and the membrane was blotted with HN1 antibody diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



An ELISA plate is coated with 50 μ L of HN1 recombinant protein at concentration ranged from 0.004 μ g/mL to 1 μ g/mL. The coated protein is detected with HN1 antibody (NBP1-32510) at 0.25 μ g/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-32510

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-32510

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

