

Product Datasheet

PKC gamma Antibody - BSA Free NBP1-32535

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-32535

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-32535



NBP1-32535

PKC gamma Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	0.1M Tris (pH 7), 0.1M Glycine, 10% Glycerol
Target Molecular Weight	78 kDa

Product Description


Host	Rabbit
Gene ID	5582
Gene Symbol	PRKCG
Species	Human, Mouse, Rat
Immunogen	Recombinant protein encompassing a sequence within the center region of human PKC gamma. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000, Immunohistochemistry-Frozen 1:100-1:1000

Images


Western Blot: PKC gamma Antibody [NBP1-32535] - Various tissue extracts (50 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with PKC gamma antibody diluted at 1:3000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.




Immunohistochemistry-Frozen: PKC gamma Antibody [NBP1-32535] - Frozen-sectioned adult mouse cerebellum. Green: PKC gamma protein stained by PKC gamma antibody diluted at 1:250. Red: beta Tubulin 3/ TUJ1, stained by beta Tubulin 3/ TUJ1 antibody [11710] diluted at 1:500. Blue: Fluoroshield with DAPI.




Immunohistochemistry-Paraffin: PKC gamma Antibody [NBP1-32535] - Paraffin-embedded Hun7 xenograft, using antibody at 1:500 dilution.




Immunohistochemistry-Paraffin: PKC gamma Antibody [NBP1-32535] - Paraffin-embedded rat brain. PKC gamma antibody diluted at 1:500.




Immunohistochemistry-Paraffin: PKC gamma Antibody [NBP1-32535] - Paraffin-embedded mouse brain. PKC gamma antibody diluted at 1:500.




Western Blot: PKC gamma Antibody [NBP1-32535] - Various tissue extracts (50 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with PKC gamma antibody (NBP1-32535) diluted at 1:3000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



Western Blot: PKC gamma Antibody [NBP1-32535] - Human tissue extract (30 ug) was separated by 7.5% SDS-PAGE, and the membrane was blotted with PKC gamma antibody (NBP1-32535) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



Western Blot: PKC gamma Antibody [NBP1-32535] - Whole cell extract (30 ug) was separated by 7.5% SDS-PAGE, and the membrane was blotted with PKC gamma antibody (NBP1-32535) diluted at 1:2000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-32535

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-38728PEP	PKC gamma Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-32535

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

