

Product Datasheet

GSTT1 Antibody - BSA Free

NBP1-32733

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-32733

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-32733



NBP1-32733

GSTT1 Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol
Target Molecular Weight	27 kDa

Product Description


Host	Rabbit
Gene ID	2952
Gene Symbol	GSTT1
Species	Human, Mouse, Rat
Reactivity Notes	Bovine (86%).
Immunogen	Full length human GSTT1 Recombinant protein.

Product Application Details


Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunoprecipitation 1:100-1:500, Immunohistochemistry-Paraffin 1:100-1:1000

Images


Immunohistochemistry-Paraffin: GSTT1 Antibody [NBP1-32733] - Human hepatoma. GSTT1 antibody dilution: 1:500. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.




Western Blot: GSTT1 Antibody [NBP1-32733] - A. 50 ug rat liver lysate/extract 12 % SDS-PAGE GSTT1 antibody dilution: 1:500




Western Blot: GSTT1 Antibody [NBP1-32733] - GSTT1 antibody detects GSTT1 protein by Western blot analysis. Whole cell extracts (30 ug) was separated by 12% SDS-PAGE, and the membrane was blotted with GSTT1 antibody (NBP1-32733) diluted by 1:1000.



Western Blot: GSTT1 Antibody [NBP1-32733] - GSTT1 antibody detects GSTT1 protein by Western blot analysis. Mouse tissue extracts (50 ug) was separated by 12 % SDS-PAGE, and the membrane was blotted with GSTT1 antibody (NBP1-32733) diluted by 1:1000.



Immunoprecipitation: GSTT1 Antibody [NBP1-32733] - Immunoprecipitation of GSTT1 protein from HepG2 whole cell extracts using 5 ug of GSTT1 antibody. Western blot analysis was performed using GSTT1 antibody. EasyBlot anti-Rabbit IgG was used as a secondary reagent.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-32733

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBC1-28782-0.1mg	Recombinant Human GSTT1 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-32733

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

