

Product Datasheet

MCM4 Antibody - BSA Free

NBP1-33105

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-33105

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-33105



NBP1-33105

MCM4 Antibody - BSA Free


Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	0.1M Tris (pH 7), 0.1M Glycine, 10% Glycerol

Product Description	
Host	Rabbit
Gene ID	4173
Gene Symbol	MCM4
Species	Human, Mouse
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human MCM4. The exact sequence is proprietary.


Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP)
Recommended Dilutions	Western Blot 1:500-1:10000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunoprecipitation 1:100-1:500, Immunohistochemistry-Paraffin 1:100-1:1000, Chromatin Immunoprecipitation (ChIP)

Images


Western Blot: MCM4 Antibody [NBP1-33105] - Various whole cell extracts (30 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with MCM4 antibody diluted at 1:1000.




Immunocytochemistry/Immunofluorescence: MCM4 Antibody [NBP1-33105] - Paraformaldehyde-fixed HeLa, using antibody at 1:500 dilution.




Immunohistochemistry-Paraffin: MCM4 Antibody [NBP1-33105] - Paraffin-embedded A549 Xenograft, using antibody at 1:100 dilution.




Western Blot: MCM4 Antibody [NBP1-33105] - MCM4 antibody detects MCM4 protein by Western blot analysis. A. 30 ug Jurkat whole cell lysate/extract. B. 30 ug Raji whole cell lysate/extract. C. 30 ug NCI-H929 whole cell lysate/extract. 7.5 % SDS-PAGE. MCM4 antibody dilution: 1:5000.




Immunoprecipitation: MCM4 Antibody [NBP1-33105] - MCM4 antibody immunoprecipitates MCM4 protein in IP experiments. IP samples: Jurkat whole cell extract. A. 40 ug Jurkat whole cell extract. B. Control with 4 ug of preimmune Rabbit IgG. C. Immunoprecipitation of MCM4 protein by 4 ug MCM4 antibody. 7.5 % SDS-PAGE. The immunoprecipitated MCM4 protein was detected by MCM4 antibody diluted at 1:500.



Immunohistochemistry: MCM4 Antibody [NBP1-33105] - Immunohistochemical analysis of paraffin-embedded zebrafish tissue, using MCM4 antibody (NBP1-33105) at 1:300 dilution.



Chromatin Immunoprecipitation (ChIP): MCM4 Antibody [NBP1-33105] - Cross-linked ChIP was performed with HeLa chromatin extract and 5 ug of either control rabbit IgG or anti-MCM4 antibody. The precipitated DNA was detected by PCR with primer set targeting to Myc promoter.




Western Blot: MCM4 Antibody [NBP1-33105] - Sample (30 ug of whole cell lysate)

A: zebrafish eye

7.5% SDS PAGE

NBP1-33105 diluted at 1:1000





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-33105

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00004173-P01-10ug	Recombinant Human MCM4 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-33105

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

