

# Product Datasheet

## TXNDC3 Antibody - BSA Free

### NBP1-33142

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-33142](http://www.novusbio.com/NBP1-33142)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-33142](http://www.novusbio.com/reviews/destination/NBP1-33142)



**NBP1-33142**

TXNDC3 Antibody - BSA Free

**Product Information**

|                                |                                                                                                                  |
|--------------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>               | 100 ul                                                                                                           |
| <b>Concentration</b>           | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>                 | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                     |
| <b>Clonality</b>               | Polyclonal                                                                                                       |
| <b>Preservative</b>            | 0.01% Thimerosal                                                                                                 |
| <b>Isotype</b>                 | IgG                                                                                                              |
| <b>Purity</b>                  | Antigen Affinity-purified                                                                                        |
| <b>Buffer</b>                  | 0.1M Tris (pH 7), 0.1M Glycine, 10% Glycerol                                                                     |
| <b>Target Molecular Weight</b> | 67 kDa                                                                                                           |

**Product Description**

|                    |                                                                                                                                                                                                |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Novus Biologicals Rabbit TXNDC3 Antibody - BSA Free (NBP1-33142) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| <b>Host</b>        | Rabbit                                                                                                                                                                                         |
| <b>Gene ID</b>     | 51314                                                                                                                                                                                          |
| <b>Gene Symbol</b> | NME8                                                                                                                                                                                           |
| <b>Species</b>     | Human                                                                                                                                                                                          |
| <b>Immunogen</b>   | Recombinant protein encompassing a sequence within the center region of human TXNDC3. The exact sequence is proprietary.                                                                       |


**Product Application Details**

|                              |                                                                                                          |
|------------------------------|----------------------------------------------------------------------------------------------------------|
| <b>Applications</b>          | Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry                                        |
| <b>Recommended Dilutions</b> | Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000 |




## Images

Western Blot: TXNDC3 Antibody [NBP1-33142] - Sample (30 ug of whole cell lysate) A: Raji 7.5% SDS PAGE; antibody diluted at 1:1000.



Immunohistochemistry-Paraffin: TXNDC3 Antibody [NBP1-33142] - Paraffin-embedded Cal27 xenograft. TXNDC3 antibody [N3C3] dilution: 1:500.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-33142**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00051314-P01-10ug | Recombinant Human TXNDC3 GST (N-Term) Protein       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-33142](http://www.novusbio.com/reviews/submit/NBP1-33142)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

