

# Product Datasheet

## CKS2 Antibody - BSA Free NBP1-33404

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-33404](http://www.novusbio.com/NBP1-33404)

Updated 9/25/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-33404](http://www.novusbio.com/reviews/destination/NBP1-33404)



**NBP1-33404**


CKS2 Antibody - BSA Free

| Product Information         |                                                                                                                                                                                                 |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 100 ul                                                                                                                                                                                          |
| Concentration               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.                                                                                |
| Storage                     | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                                                                                                    |
| Clonality                   | Polyclonal                                                                                                                                                                                      |
| Preservative                | 0.01% Thimerosal                                                                                                                                                                                |
| Isotype                     | IgG                                                                                                                                                                                             |
| Purity                      | Antigen Affinity-purified                                                                                                                                                                       |
| Buffer                      | 0.1M Tris, 0.1M Glycine, 10% Glycerol                                                                                                                                                           |
| Target Molecular Weight     | 10 kDa                                                                                                                                                                                          |
| Product Description         |                                                                                                                                                                                                 |
| Description                 | Novus Biologicals Rabbit CKS2 Antibody - BSA Free (NBP1-33404) is a polyclonal antibody validated for use in WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host                        | Rabbit                                                                                                                                                                                          |
| Gene ID                     | 1164                                                                                                                                                                                            |
| Gene Symbol                 | CKS2                                                                                                                                                                                            |
| Species                     | Human                                                                                                                                                                                           |
| Reactivity Notes            | Xenopus laevis (100%).                                                                                                                                                                          |
| Immunogen                   | Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human CKS2. The exact sequence is proprietary.                                             |
| Product Application Details |                                                                                                                                                                                                 |
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence                                                                                                                                           |
| Recommended Dilutions       | Western Blot 1:500-1:3000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000                                                                                                                 |




## Images

Immunocytochemistry/Immunofluorescence: CKS2 Antibody [NBP1-33404] - Paraformaldehyde-fixed HeLa, using antibody at 1:500 dilution.



Whole cell extract (30 ug) was separated by 15% SDS-PAGE, and the membrane was blotted with CKS2 antibody [N2C3] (NBP1-33404) diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP1-33404**

---

|             |                                                     |
|-------------|-----------------------------------------------------|
| NBP2-33376H | Blue Marker Antibody (6F4-F6) [HRP]                 |
| HAF008      | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160      | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891  | Rabbit IgG Isotype Control                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-33404](http://www.novusbio.com/reviews/submit/NBP1-33404)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

