

Product Datasheet

ATP Citrate Lyase Antibody - BSA Free NBP1-33430

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-33430

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-33430




NBP1-33430**ATP Citrate Lyase Antibody - BSA Free**

Product Information	
Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	0.1M Tris (pH 7), 0.1M Glycine, 10% Glycerol
Target Molecular Weight	121 kDa
Product Description	
Host	Rabbit
Gene ID	47
Gene Symbol	ACLY
Species	Human, Mouse, Rat
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human ATP Citrate Lyase. The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:500-1:3000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000




Images


Western Blot: ATP Citrate Lyase Antibody [NBP1-33430] - Sample (30 ug of whole cell lysate) A: HeLa B: Raji 7. 5% SDS PAGE; antibody diluted at 1:1000.




Immunocytochemistry/Immunofluorescence: ATP Citrate Lyase Antibody [NBP1-33430] - HeLa cells were fixed in ice-cold MeOH for 5 min. Green: ATP citrate lyase protein stained by ATP citrate lyase antibody [C2C3], C-term diluted at 1:200. Blue: Hoechst 33342 staining.



Western Blot: ATP Citrate Lyase Antibody [NBP1-33430] - Sample (50 ug of whole cell lysate) A: mouse brain 5% SDS PAGE, antibody diluted at 1:2000.



Western Blot: ATP Citrate Lyase Antibody [NBP1-33430] - A. 30 ug PC-12 whole cell lysate/extract, B. 30 ug Rat 2 whole cell lysate/extract 5% SDS-PAGE gel, antibody dilution 1:1000.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-33430

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-90269PEP	ATP Citrate Lyase Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-33430

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

