

# Product Datasheet

## PD-1 Antibody (J116) [Allophycocyanin] NBP1-43107APC

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-43107APC](http://www.novusbio.com/NBP1-43107APC)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-43107APC](http://www.novusbio.com/reviews/destination/NBP1-43107APC)



**NBP1-43107APC**

PD-1 Antibody (J116) [Allophycocyanin]

| Product Information |  |
|---------------------|--|
| Unit Size           | 0.1 ml   |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.   |
| Clonality           | Monoclonal   |
| Clone               | J116   |
| Preservative        | 0.05% Sodium Azide   |
| Isotype             | IgG1 Kappa   |
| Conjugate           | Allophycocyanin  |
| Purity              | Protein A or G purified  |
| Buffer              | PBS  |

| Product Description |  |
|---------------------|--|
| Host                | Mouse                                    |
| Gene ID             | 5133                                     |
| Gene Symbol         | PDCD1                                    |
| Species             | Human                                    |
| Immunogen           | The immunogen for this antibody was PD1. |

| Product Application Details |  |
|-----------------------------|--|
| Applications                | Western Blot, Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunohistochemistry, Immunoprecipitation, Immunohistochemistry-Frozen |

**Images**

Product Image: PD-1 Antibody (J116) [Allophycocyanin] [NBP1-43107APC] - Vial of APC conjugated antibody. APC is optimally excited at 650 nm by the Red laser (633 or 640 nm) and has an emission maximum of 660 nm.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-43107APC**

---

|               |  |
|---------------|--|
| IC002A        | Mouse IgG1 Isotype Control (11711) [Allophycocyanin] |
| NBP1-43107V   | PD-1 Antibody (J116) [DyLight 405]                   |
| NBP1-77276PEP | PD-1 Antibody Blocking Peptide                       |
| 210-TA-005    | TNF-alpha [Unconjugated]                             |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-43107APC](http://www.novusbio.com/reviews/submit/NBP1-43107APC)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

