

Product Datasheet

Aminopeptidase A/ENPEP Antibody (6C3) - BSA Free NBP1-43281-0.5mg

Unit Size: 0.5 mg

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-43281

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-43281



NBP1-43281-0.5mg

Aminopeptidase A/ENPEP Antibody (6C3) - BSA Free

Product Information	
Unit Size	0.5 mg
Concentration	1.0 mg/ml
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	6C3
Preservative	0.09% Sodium Azide
Isotype	IgG2a Kappa
Purity	Protein A or G purified
Buffer	PBS


Product Description	
Host	Rat
Gene ID	2028
Gene Symbol	ENPEP
Species	Mouse
Immunogen	The immunogen for this antibody was Aminopeptidase A/APA.

Product Application Details	
Applications	Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Flow Cytometry 1:10-1:1000, Immunohistochemistry 1:10-1:500, Immunoprecipitation 1:10-1:500, Immunohistochemistry-Paraffin 1:10-1:500




Images

Flow Cytometry: Aminopeptidase A/ENPEP Antibody (6C3) [NBP1-43281] - Detection of Aminopeptidase A/ENPEP in Mouse Bone Marrow Cells by Flow Cytometry. Mouse bone marrow cells were stained with PE-conjugated Rat Anti-Mouse CD11b (Catalog # FAB1124P) and Rat Anti-Mouse Aminopeptidase A/ENPEP Monoclonal Antibody (Catalog # NBP1-43281, top panel) or isotype control antibody (Catalog # MAB006, bottom panel), followed by APC-conjugated Goat Anti-Rat IgG (Catalog # F0113).



Flow Cytometry: Aminopeptidase A/ENPEP Antibody (6C3) [NBP1-43281] - Staining of BALB/c bone marrow cells with 0.5 ug of Rat IgG2a K Isotype Control Purified (open histogram) or 0.5 ug of Anti-Mouse BP-1 Purified (filled histogram) followed by Anti-Rat IgG Biotin and Streptavidin PE . Cells in the lymphocyte gate were used





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-43281-0.5mg

HAF005	Goat anti-Rat IgG Secondary Antibody [HRP]
NBP1-75398	Goat anti-Rat IgG (H+L) Secondary Antibody (Pre-adsorbed)
NBP1-43321-0.5mg	Rat IgG2a Kappa Light Chain Isotype Control (R2a)
NBP1-43281AF405	Aminopeptidase A/ENPEP Antibody (6C3) [Alexa Fluor® 405]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-43281

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

