

Product Datasheet

PDGFR alpha Antibody (APA5) [DyLight 650] NBP1-43350C

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-43350C

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-43350C



NBP1-43350C

PDGFR alpha Antibody (APA5) [DyLight 650]

| Product Information | |
|------------------------------------|--|
| Unit Size | 0.1 ml |
| Concentration | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | APA5 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2a Kappa |
| Conjugate | DyLight 650 |
| Purity | Protein A or G purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Rat |
| Gene ID | 5156 |
| Gene Symbol | PDGFRA |
| Species | Human, Mouse |
| Reactivity Notes | Human reactivity reported from a verified customer review. |
| Immunogen | The immunogen for this antibody was CD140a. |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |
| Product Application Details | |
| Applications | Western Blot, Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen |
| Recommended Dilutions | Western Blot, Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. In flow cytometry, we recommend using a cell surface staining protocol and fixation and permeabilization is not required to get staining. If fixation/permeabilization is needed for other markers in the same experiment, we recommend doing all surface staining prior to the fix/perm procedure. This antibody may not bind to a fixed epitope. |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-43350C

| | |
|--------------------|---|
| NBP1-51104C-0.25ml | Rat IgG2a Isotype Control (KLH/G2a-1-1) [DyLight 650] |
| NBP2-76508PEP | PDGFR alpha Recombinant Protein Antigen |
| DVE00 | VEGF [HRP] |
| DYC2114-2 | PDGFR alpha [Biotin] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-43350C

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

