

Product Datasheet

CD117/c-kit Antibody (ACK2) [Allophycocyanin] NBP1-43359APC

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-43359APC

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-43359APC



NBP1-43359APC

CD117/c-kit Antibody (ACK2) [Allophycocyanin]

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	ACK2
Preservative	0.05% Sodium Azide
Isotype	IgG2b Kappa
Conjugate	Allophycocyanin
Purity	Protein A or G purified
Buffer	PBS
Product Description	
Host	Rat
Gene ID	3815
Gene Symbol	KIT
Species	Human, Mouse, Guinea Pig
Reactivity Notes	Human reactivity reported in scientific literature (PMID: 23967247). Guinea Pig reactivity reported in scientific literature (PMID: 21654812)
Marker	Hematopoietic Stem Cell Marker
Specificity/Sensitivity	The ACK2 monoclonal antibody reacts with mouse CD117, also known as c-Kit receptor, Steel factor receptor and stem cell factor receptor.
Immunogen	Mouse CD117
Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation, Immunohistochemistry Whole-Mount
Recommended Dilutions	Western Blot, Flow Cytometry, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Immunohistochemistry-Frozen, Immunohistochemistry Whole-Mount
Application Notes	Optimal dilution of this antibody should be experimentally determined.



Images

CD117/c-kit Antibody (ACK2) [Allophycocyanin] [NBP1-43359APC] - Vial of APC conjugated antibody. APC is optimally excited at 650 nm by the Red laser (633 or 640 nm) and has an emission maximum of 660 nm.



Allophycocyanin (APC)

LASER (nm)	FILTER
Red (633,640)	660/10

EXCITATION MAX (nm)	EMISSION MAX (nm)
650	660



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-43359APC

NBP1-43323APC	Rat IgG2b Kappa Light Chain Isotype Control (149/10H5) [Allophycocyanin]
NBP2-52975PEP	CD117/c-kit Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]
DY332	CD117/c-kit [Biotin]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-43359APC

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

