

# Product Datasheet

## Integrin alpha E/CD103 Antibody (B-Ly7) [DyLight 405] NBP1-43370V

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-43370V](http://www.novusbio.com/NBP1-43370V)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-43370V](http://www.novusbio.com/reviews/destination/NBP1-43370V)



**NBP1-43370V**

Integrin alpha E/CD103 Antibody (B-Ly7) [DyLight 405]

| Product Information         |  |
|-----------------------------|--|
| Unit Size                   | 0.1 ml   |
| Concentration               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage                     | Store at 4C in the dark.   |
| Clonality                   | Monoclonal   |
| Clone                       | B-Ly7  |
| Preservative                | 0.05% Sodium Azide   |
| Isotype                     | IgG1 Kappa   |
| Conjugate                   | DyLight 405  |
| Purity                      | Protein A or G purified  |
| Buffer                      | 50mM Sodium Borate   |
| Product Description         |  |
| Host                        | Mouse  |
| Gene ID                     | 3682   |
| Gene Symbol                 | ITGAE  |
| Species                     | Human  |
| Immunogen                   | The immunogen for this antibody was CD103/Integrin alpha E.  |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                |
| Product Application Details |  |
| Applications                | Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen  |
| Recommended Dilutions       | Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.   |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-43370V**

|                    |   |
|--------------------|---|
| NBP1-43319V-0.5ml  | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1) [DyLight 405]   |
| H00003682-Q01-10ug | Recombinant Human Integrin alpha E/CD103 GST (N-Term) Protein |
| 210-TA-005         | TNF-alpha [Unconjugated]                                      |
| NBP2-10814         | Integrin alpha E/CD103 Overexpression Lysate                  |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-43370V](http://www.novusbio.com/reviews/submit/NBP1-43370V)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

