

Product Datasheet

B7-1/CD80 Antibody (16-10A1) [Allophycocyanin/Cy7] NBP1-43385APCCY7

Unit Size: 0.1 ml

Store at 4C in the dark. Do not freeze.

www.novusbio.com



technical@novusbio.com


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-43385APCCY7

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-43385APCCY7



NBP1-43385APCCY7

B7-1/CD80 Antibody (16-10A1) [Allophycocyanin/Cy7]

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark. Do not freeze.
Clonality	Monoclonal
Clone	16-10A1
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	Allophycocyanin/Cy7
Purity	Protein A or G purified
Buffer	PBS

Product Description	
Host	Armenian Hamster
Gene ID	941
Gene Symbol	CD80
Species	Human, Mouse, Feline
Reactivity Notes	Feline reactivity reported from a verified customer review.
Immunogen	The immunogen for this antibody was CD80.

Product Application Details	
Applications	Flow Cytometry, Product Image
Recommended Dilutions	Flow Cytometry, Product Image
Application Notes	Optimal dilution of this antibody should be experimentally determined. For optimal results using our Tandem dyes, please avoid prolonged exposure to light or extreme temperature fluctuations. These can lead to irreversible degradation or decoupling. When staining intracellular targets, specific attention to the fixation and permeabilization steps in your flow protocol may be required. Please contact our technical support team at technical@novusbio.com if you have any questions.

Images

Product Image: B7-1/CD80 Antibody (16-10A1) [Allophycocyanin/Cy7] [NBP1-43385APCCY7] - Vial of APC/Cy7 conjugated antibody. APC/Cy7 is optimally excited at 652 nm by the Red laser (633 or 640 nm) and has an emission maximum of 775 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-43385APCCY7

NBP2-26578	Recombinant Mouse B7-1/CD80 Protein
210-TA-005	TNF-alpha [Unconjugated]
740-B1-100	B7-1/CD80 [Unconjugated]
6507-IL-010/CF	IL-4 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-43385APCCY7

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

