

Product Datasheet

CD1d Antibody (51.1) [mFluor Violet 450 SE] NBP1-43460MFV450

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-43460MFV450

Updated 9/20/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-43460MFV450



NBP1-43460MFV450

CD1d Antibody (51.1) [mFluor Violet 450 SE]

| Product Information | |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size | 0.1 ml |
| Concentration | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 51.1 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2b |
| Conjugate | mFluor Violet 450 SE |
| Purity | Protein A or G purified |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|---------------------|-------------------------------------------|
| Host | Mouse |
| Gene ID | 912 |
| Gene Symbol | CD1D |
| Species | Human |
| Immunogen | The immunogen for this antibody was CD1d. |

| Product Application Details | |
|-----------------------------|----------------------------------------------------------------------------------------|
| Applications | Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation |
| Recommended Dilutions | Flow Cytometry, Immunohistochemistry, Immunoprecipitation, Immunohistochemistry-Frozen |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |

Images

CD1d Antibody (51.1) [mFluor Violet 450 SE] [NBP1-43460MFV450] - Vial of mFluor Violet 450 conjugated antibody. mFluor Violet 450 is optimally excited at 406 nm by the Violet laser (405 nm) and has an emission maximum of 445 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-43460MFV450

| | |
|-----------------|------------------------------------|
| NBP2-61190-10ug | Recombinant Mouse CD1d His Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| NBL1-08914 | CD1d Overexpression Lysate |
| 6507-IL-010/CF | IL-4 [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-43460MFV450

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

