

Product Datasheet

USP13 Antibody (OTI4G11)

NBP1-48047

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-48047

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-48047



NBP1-48047

USP13 Antibody (OTI4G11)

Product Information	
Unit Size	0.1 ml
Concentration	0.41 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	OTI4G11
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.3), 1.0% BSA and 50% Glycerol
Target Molecular Weight	97.3 kDa


Product Description	
Description	Novus Biologicals Mouse USP13 Antibody (OTI4G11) (NBP1-48047) is a monoclonal antibody validated for use in IHC, WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	8975
Gene Symbol	USP13
Species	Human, Mouse, Rat, Canine, Monkey, Primate
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Mouse-On-Mouse blocking reagent may be needed for IHC and ICC experiments to reduce high background signal. You can find these reagents under catalog numbers PK-2200-NB and MP-2400-NB. Please contact Technical Support if you have any questions.
Specificity/Sensitivity	This antibody is specific for Homo sapiens ubiquitin specific peptidase 13 (isopeptidase T-3) (USP13).
Immunogen	Full-length protein expressed in 293T cell transfected with human USP13 expression vector

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000, Immunohistochemistry 1:50, Immunoprecipitation 4ug/ml, Immunohistochemistry-Paraffin 1:50




Images


Western Blot: USP13 Antibody (OTI4G11) [NBP1-48047] - Analysis of extracts (35ug) from 9 different cell lines by using anti-USP13 monoclonal antibody.




Immunohistochemistry-Paraffin: USP13 Antibody (OTI4G11) [NBP1-48047] - Staining of paraffin-embedded prostate tissue using anti-USP13 mouse monoclonal antibody.




Western Blot: USP13 Antibody (OTI4G11) [NBP1-48047] - HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY USP13 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-USP13.




Immunohistochemistry-Paraffin: USP13 Antibody (OTI4G11) [NBP1-48047] - Staining of paraffin-embedded breast tissue using anti-USP13 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: USP13 Antibody (OTI4G11) [NBP1-48047] - Staining of paraffin-embedded Carcinoma of liver tissue using anti-USP13 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: USP13 Antibody (OTI4G11) [NBP1-48047] - Staining of paraffin-embedded colon tissue using anti-USP13 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: USP13 Antibody (OTI4G11) [NBP1-48047] - Staining of paraffin-embedded Kidney tissue using anti-USP13 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: USP13 Antibody (OTI4G11) [NBP1-48047] - Staining of paraffin-embedded liver tissue using anti-USP13 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: USP13 Antibody (OTI4G11) [NBP1-48047] - Staining of paraffin-embedded lung tissue using anti-USP13 mouse monoclonal antibody.




Immunohistochemistry-Paraffin: USP13 Antibody (OTI4G11) [NBP1-48047] - Staining of paraffin-embedded Ovary tissue using anti-USP13 mouse monoclonal antibody.



Immunohistochemistry-Paraffin: USP13 Antibody (OTI4G11) [NBP1-48047] - Staining of paraffin-embedded pancreas tissue using anti-USP13 mouse monoclonal antibody.



Immunoprecipitation: USP13 Antibody (OTI4G11) [NBP1-48047] - (Negative control: IP without adding anti-USP13 antibody.). For each experiment, 500ul of DDK tagged USP13 overexpression lysates (at 1:5 dilution with HEK293T lysate), 2ug of anti-USP13 antibody and 20ul (0.1mg) of goat anti-mouse conjugated magnetic beads were mixed and incubated overnight. After extensive wash to remove any non-specific binding, the immuno-precipitated products were analyzed with rabbit anti-DDK polyclonal antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-48047

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP1-83088PEP	USP13 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-48047

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

