

Product Datasheet

JAM-A Antibody NBP1-49908

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-49908

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-49908



NBP1-49908


JAM-A Antibody

Product Information	
Unit Size	0.1 ml
Concentration	0.2 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	TBS and 0.1% BSA
Product Description	
Host	Rabbit
Gene ID	50848
Gene Symbol	F11R
Species	Human, Mouse
Reactivity Notes	Use in Mouse reported in scientific literature (PMID:34293456)
Immunogen	The immunogen for this product maps to a region between residue 249 and 299 of human Junctional Adhesion Molecule A using the numbering given in entry NP_058642.1 (GeneID 50848).
Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Western Blot 1:2000-1:10000, Immunohistochemistry 1:100-1:500, Immunoprecipitation 2-5 ug/mg lysate, Immunohistochemistry-Paraffin 1:100-1:500
Application Notes	Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.




Images

Western Blot: JAM-A Antibody [NBP1-49908] - Samples: Whole cell lysate from 293T (5, 15 and 50 ug for WB; 1 mg for IP, 20% of IP loaded) and HeLa (H; 50 ug) cells. Antibodies: Affinity purified rabbit anti-JAM-A antibody NBP1-49908 used for WB at 0.04 ug/ml (A) and 1 ug/ml (B) and used for IP at 3 ug/mg lysate. JAM-A was not successfully immunoprecipitated by rabbit anti-JAM-A antibody NBP1-49907, which recognizes an upstream epitope. Detection: Chemiluminescence with exposure times of 3 minutes (A and B).



Immunohistochemistry-Paraffin: JAM-A Antibody [NBP1-49908] - Human prostate carcinoma. Antibody: Affinity purified rabbit anti-JAM-A used at a dilution of 1:200 (1ug/ml). Detection: DAB



Publications

Pan Z, Gao Y, Liu S et al. Wu-Zi-Yan-Zong-Wan protects mouse blood-testis barrier from Tripterygium wilfordii Hook. f. multiglycoside-induced disruption by regulating proinflammatory cytokines Journal of ethnopharmacology 2021-07-19 [PMID: 34293456] (WB)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-49908

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-38934PEP	JAM-A Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-49908

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

