

Product Datasheet

Casein Kinase 1 delta Antibody - BSA Free NBP1-52975

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-52975

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-52975



NBP1-52975

Casein Kinase 1 delta Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS, 2% Sucrose
Target Molecular Weight	47 kDa

Product Description

Description	The addition of 50% glycerol is optional for those storing this antibody at -20C and not aliquoting smaller units. However, please note that glycerol may interrupt some downstream antibody applications and should be added with caution.
Host	Rabbit
Gene ID	1453
Gene Symbol	CSNK1D
Species	Human
Immunogen	Synthetic peptides corresponding to CSNK1D(casein kinase 1, delta) The peptide sequence was selected from the middle region of CSNK1D. Peptide sequence NLVYIIDFGLAKKYRDARTHQHIPPYRENKNLTGTARYASINTHLGIEQSR. The peptide sequence for this immunogen was taken from within the described region.

Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot 1.0 ug/ml



Images

Western Blot: Casein Kinase 1 delta Antibody [NBP1-52975] - Sample Type: Human Adult Placenta Antibody Dilution: 1.0 ug/ml

CSNK1D

90—
65—
40—
31—
22—

Rabbit Anti-CSNK1D
Sample Type: Human Adult Placenta
Antibody Concentration: 1ug/mL

Western Blot: Casein Kinase 1 delta Antibody [NBP1-52975] - OVCAR-3 cell lysate, concentration 0.2-1 ug/ml.

90 kDa—
65 kDa—
40 kDa—
31 kDa—
22 kDa—

Western Blot: Casein Kinase 1 delta Antibody [NBP1-52975] - Human Fetal Lung, Antibody Dilution: 1.0 ug/ml.

90—
65—
40—
31—
22—



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-52975

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-51703-0.1mg	Recombinant Human Casein Kinase 1 delta His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-52975

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

