

# Product Datasheet

## FKBP25 Antibody - BSA Free

### NBP1-54354

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-54354](http://www.novusbio.com/NBP1-54354)

Updated 2/21/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-54354](http://www.novusbio.com/reviews/destination/NBP1-54354)



**NBP1-54354**

FKBP25 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS, 2% Sucrose
Target Molecular Weight	25 kDa
Product Description	
Description	The addition of 50% glycerol is optional for those storing this antibody at -20C and not aliquoting smaller units. However, please note that glycerol may interrupt some downstream antibody applications and should be added with caution.
Host	Rabbit
Gene ID	2287
Gene Symbol	FKBP3
Species	Human
Immunogen	Synthetic peptides corresponding to FKBP3(FK506 binding protein 3, 25kDa) The peptide sequence was selected from the C terminal of FKBP3. Peptide sequence EALLTMSKGEKARLEIEPEWAYGKKKGQPDAKIPPNAKLTFEVELVDID. The peptide sequence for this immunogen was taken from within the described region.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1.0 ug/ml



## Images


Western Blot: FKBP25 Antibody [NBP1-54354] - Human Hela, Antibody Dilution: 1.0 ug/ml FKBP3 is supported by BioGPS gene expression data to be expressed in HeLa.




Western Blot: FKBP25 Antibody [NBP1-54354] - Jurkat cell lysate, concentration 0.2-1 ug/ml.



Western Blot: FKBP25 Antibody [NBP1-54354] - Human 721\_B, Antibody Dilution: 1.0 ug/ml FKBP3 is supported by BioGPS gene expression data to be expressed in 721\_B.



Western Blot: FKBP25 Antibody [NBP1-54354] - Human MCF7, Antibody Dilution: 1.0 ug/ml FKBP3 is supported by BioGPS gene expression data to be expressed in MCF7.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-54354**

---

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-41160	Recombinant Human FKBP25 Protein

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-54354](http://www.novusbio.com/reviews/submit/NBP1-54354)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

