

Product Datasheet

WDR33 Antibody - BSA Free

NBP1-60030

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-60030

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-60030



NBP1-60030

WDR33 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS, 2% Sucrose

Product Description	
Description	The addition of 50% glycerol is optional for those storing this antibody at -20C and not aliquoting smaller units. However, please note that glycerol may interrupt some downstream antibody applications and should be added with caution.
Host	Rabbit
Gene ID	55339
Gene Symbol	WDR33
Species	Human
Immunogen	Synthetic peptides corresponding to WDR33(WD repeat domain 33) The peptide sequence was selected from the middle region of WDR33. Peptide sequence TKFVRTSTNKVKCPVFVVRWTPEGRRRLVTGASSGEFTLWNGLTFFNFETIL. The peptide sequence for this immunogen was taken from within the described region.


Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1.0 ug/ml



Images

Western Blot: WDR33 Antibody [NBP1-60030] - Human HepG2, Antibody Dilution: 1.0 ug/ml WDR33 is supported by BioGPS gene expression data to be expressed in HepG2.

WDR33




Western Blot: WDR33 Antibody [NBP1-60030] - HepG2 cell lysate, concentration 0.2-1 ug/ml.



Western Blot: WDR33 Antibody [NBP1-60030] - Lanes: 1 : SREC pulldown from lysate from 10^6 human 293T cells Primary, Antibody Dilution: 1 : 1000 Secondary Antibody: Anti-rabbit-Alexa Fluor Secondary, Antibody Dilution: 1 : 5000 Gene name: WDR33.

WDR33



See Immunoblot 2 Data and customer Feedback for more Information



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-60030

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-83756PEP	WDR33 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-60030

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

