

# Product Datasheet

## CD31/PECAM-1 Antibody - BSA Free NBP1-71663

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

**Reviews: 2 Publications: 17**

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-71663](http://www.novusbio.com/NBP1-71663)

Updated 2/17/2025 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-71663](http://www.novusbio.com/reviews/destination/NBP1-71663)




**NBP1-71663**

CD31/PECAM-1 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS, 30% Glycerol
Target Molecular Weight	82.5 kDa
Product Description	
Host	Rabbit
Gene ID	5175
Gene Symbol	PECAM1
Species	Human, Mouse
Immunogen	This CD31/PECAM-1 Antibody was developed against a genomic peptide made to an internal region of human CD31/PECAM1 (within residues 100-300) [Swiss-Prot P16284].
Notes	Manufactured by Genomic Antibody Technology™. GAT <a href="#">FAQs</a>
Product Application Details	
Applications	Western Blot, Simple Western, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:10000, Simple Western 1:1000, Flow Cytometry 1:1000, Immunohistochemistry 1:50, Immunocytochemistry/ Immunofluorescence 1:10-1:500, Immunohistochemistry-Paraffin 1:50, Immunohistochemistry-Frozen reported in scientific literature (PMID 35879332)
Application Notes	<p>In Western Blot a band is seen ~150 kDa. Prior to immunostaining paraffin tissues, antigen retrieval with sodium citrate buffer (pH 6.0) is recommended. Customer feedback has been negative on cryosections.</p> <p>In Simple Western only 10 - 15 uL of the recommended dilution is used per data point. See <a href="#">Simple Western Antibody Database</a> for Simple Western validation: Tested in HUVEC lysate 0.5 mg/mL, separated by Size, antibody dilution of 1:1000, apparent MW was 166 kDa. Separated by Size-Wes, Sally Sue/Peggy Sue.</p>

## Images


Immunocytochemistry/Immunofluorescence: CD31/PECAM-1 Antibody [NBP1-71663] - Staining of mouse glioma tissue using CD31/PECAM-1 Antibody (NBP1-71663) at a dilution of 1:2500. Image provided by Dr. Eric Woolf of Dignity Health.




Simple Western: CD31/PECAM-1 Antibody [NBP1-71663] - Image shows a specific band for CD31/PECAM 1 in 0.5 mg/mL of HUVEC lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system.




Immunohistochemistry: CD31/PECAM-1 Antibody [NBP1-71663] - 5-HTR1A agonist enhances neovascularization during skin wound healing (excisional punch biopsy model). CD31/PECAM-1 endothelial cell immunohistochemical staining at day 5 (A,B) and day 7 (E,F), image scale bars represented at 50  $\mu$ m. Results were presented as mean  $\pm$  SEM (n = 5 images per wound). Two-way ANOVA was performed and significance levels were set at \*P < 0.05, \*\*P < 0.01. Image collected and cropped by CiteAb from the following publication (<https://www.frontiersin.org/article/10.3389/fphar.2018.01406/full>), licensed under a CC-BY license.




Western Blot: CD31/PECAM-1 Antibody [NBP1-71663] - Snail and EndMT-related protein expression in lungs in orthotopic murine 4T1 breast cancer model. ag Expression of Snail (see Methods). Levels of CD31/PECAM-1 determined by western blot analysis in pooled (n = 6) samples in control and 4T1 breast cancer-bearing mice from the 1st to 5th week after 4T1 cancer cell inoculation (see Methods). Results presented as fold change vs control sample corresponding to healthy mice. Total protein after transfer was used as loading control Image collected and cropped by Citeab from the following publication (Nitric oxide deficiency and endothelial-mesenchymal transition of pulmonary endothelium in the progression of 4T1 metastatic breast cancer in mice. Breast Cancer Res (2018) licensed under a CC-BY license




**Immunocytochemistry/Immunofluorescence: CD31/PECAM-1 Antibody [NBP1-71663]** - HUVEC cells were fixed for 10 minutes using 10% formalin and then permeabilized for 5 minutes using 1X TBS + 0.5% Triton X-100. The cells were incubated with CD31/PECAM-1 Antibody at 5 ug/mL overnight at 4C and detected with an anti-rabbit Dylight 488 (Green) at a 1:500 dilution. Alpha tubulin (DM1A) NB100-690 was used as a co-stain at a 1:1000 dilution and detected with an anti-mouse Dylight 550 (Red) at a 1:500 dilution. Nuclei were counterstained with DAPI (Blue). Cells were imaged using a 40X objective.



**Flow Cytometry: CD31/PECAM-1 Antibody [NBP1-71663]** - CD31/PECAM-1 Antibody was tested in WeHi cells with Dylight 488 (green) alongside a matched isotype control (black).



**Western Blot: CD31/PECAM-1 Antibody [NBP1-71663]** - Detection of CD31 in HUVEC cell lysates.



## Publications

Smeda M, Jaszta A, Maleki EH et al. Endothelial-mesenchymal transition induced by metastatic 4T1 breast cancer cells in pulmonary endothelium in aged mice *Frontiers in Molecular Biosciences* 2022-11-24 [PMID: 36504711] (In vivo assay)

Marycz K, Bourebaba N, Serwotka-Suszczak A et al. In Vitro Generated Equine Hepatic-Like Progenitor Cells as a Novel Potent Cell Pool for Equine Metabolic Syndrome (EMS) Treatment *Stem Cell Reviews and Reports* 2023-05-01 [PMID: 36658383] (Flow Cytometry, Western Blot)

Feng D, Zhou J, Liu H et al. Astrocytic NDRG2-PPM1A interaction exacerbates blood-brain barrier disruption after subarachnoid hemorrhage *Science advances* 2022-09-30 [PMID: 36179025] (ICC/IF, Mouse)

### Details:

Dilution used in ICC/IF 1:100

Bae E, Huang P, Muller-Greven G et al. Integrin alpha 3 beta 1 promotes vessel formation of glioblastoma-associated endothelial cells through calcium-mediated macropinocytosis and lysosomal exocytosis *Nature communications* 2022-07-25 [PMID: 35879332] (ICC/IF, IHC-P, IHC-Fr, Human)

Zhang L, Soni S, Hekimoglu E et al. Impaired Autophagic Activity Contributes to the Pathogenesis of Bronchopulmonary Dysplasia: Evidence from Murine and Baboon Models *Am. J. Respir. Cell Mol. Biol.* 2020-05-06 [PMID: 32374619] (WB)

Zhu NN, Lu MJ, Chen YQ et al. Autologous blood transfusion stimulates wound healing in diabetic mice through activation of the HIF-1 alpha pathway by improving the blood preservation solution *FASEB J.* 2020-03-23 [PMID: 32202355] (IHC-P, IF/IHC)

Tsai HF, Toda-Peters K, Shen AQ. Glioblastoma adhesion in a quick-fit hybrid microdevice *Biomed Microdevices* 2019-03-21 [PMID: 30900024] (ICC/IF, Human)

Sadiq A, Menchetti I, Shah A et al. 5-HT1A Receptor Function Makes Wound Healing a Happier Process *Front. Pharmacol.* 2018-12-11 [PMID: 30618734] (IHC-P, Mouse)

Smeda M, Kieronska A, et al. Nitric oxide deficiency and endothelial-mesenchymal transition of pulmonary endothelium in the progression of 4T1 metastatic breast cancer in mice. *Breast Cancer Res* 2018-08-03 [PMID: 30075800] (WB, Mouse)

Smeda M, Kieronska A, et al. Dual antiplatelet therapy with clopidogrel and aspirin increases mortality in 4T1 metastatic breast cancer-bearing mice by inducing vascular mimicry in primary tumour. *Oncotarget* 2018-04-03 [PMID: 29707148] (WB, Mouse)

Ma Y, Huang YX, Chen YY. miRNA-34a-5p downregulation of VEGFA in endometrial stem cells contributes to the pathogenesis of endometriosis. *Mol Med Rep.* 2017-09-29 [PMID: 28990049] (Human)

### Details:

This citation used the FITC form of this antibody.

Abraham V, Parambath A, Joe DS, DeLisser HM. Influence of PECAM-1 ligand interactions on PECAM-1-dependent cell motility and filopodia extension. *Physiol Rep.* 2016-11-01 [PMID: 27895229] (Human, Mouse)

### Details:

This citation used the DyLight 650 version of this antibody.

More publications at <http://www.novusbio.com/NBP1-71663>



## Procedures

### Western Blot protocolspecific for CD31 antibody (NBP1-71663)

#### Western Blot Protocol

1. Perform SDS-PAGE on samples to be analyzed, loading 40 ug of total protein per lane.
2. Transfer proteins to membrane according to the instructions provided by the manufacturer of the membrane and transfer apparatus.
3. Stain according to standard Ponceau S procedure (or similar product) to assess transfer success, and mark molecular weight standards where appropriate.
4. Rinse the blot.
5. Block the membrane using standard blocking buffer for at least 1 hour.
6. Wash the membrane in wash buffer three times for 10 minutes each.
7. Dilute primary antibody in blocking buffer and incubate 1 hour at room temperature.
8. Wash the membrane in wash buffer three times for 10 minutes each.
9. Apply the diluted HRP conjugated secondary antibody in blocking buffer (as per manufacturers instructions) and incubate 1 hour at room temperature.
10. Wash the blot in wash buffer three times for 10 minutes each (this step can be repeated as required to reduce background).
11. Apply the detection reagent of choice in accordance with the manufacturers instructions.

\*Note: Tween-20 can be added to the blocking or antibody dilution buffer at a final concentration of 0.05-0.2%.

### Immunohistochemistry-Paraffin protocol for CD31/PECAM-1 Antibody (NBP1-71663)

#### Immunohistochemistry-Paraffin Embedded Sections

##### Antigen Unmasking:

Bring slides to a boil in 10 mM sodium citrate buffer (pH 6.0) then maintain at a sub-boiling temperature for 10 minutes. Cool slides on bench-top for 30 minutes.

##### Staining:

1. Wash sections in deionized water three times for 5 minutes each.
2. Wash sections in wash buffer for 5 minutes.
3. Block each section with 100-400 ul blocking solution for 1 hour at room temperature.
4. Remove blocking solution and add 100-400 ul diluted primary antibody. Incubate overnight at 4C.
5. Remove antibody solution and wash sections in wash buffer three times for 5 minutes each.
6. Add 100-400 ul biotinylated diluted secondary antibody. Incubate 30 minutes at room temperature.
7. Remove secondary antibody solution and wash sections three times with wash buffer for 5 minutes each.
8. Add 100-400 ul Streptavidin-HRP reagent to each section and incubate for 30 minutes at room temperature.
9. Wash sections three times in wash buffer for 5 minutes each.
10. Add 100-400 ul DAB substrate to each section and monitor staining closely.
11. As soon as the sections develop, immerse slides in deionized water.
12. Counterstain sections in hematoxylin.
13. Wash sections in deionized water two times for 5 minutes each.
14. Dehydrate sections.
15. Mount coverslips.



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-71663**

---

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-71663H	CD31/PECAM-1 Antibody [HRP]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-71663](http://www.novusbio.com/reviews/submit/NBP1-71663)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

