

# Product Datasheet

## Goat anti-Mouse IgA/IgG/IgM (H+L) Secondary Antibody [DyLight 650] (Pre-adsorbed) NBP1-72561C

Unit Size: 0.5 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-72561C](http://www.novusbio.com/NBP1-72561C)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-72561C](http://www.novusbio.com/reviews/destination/NBP1-72561C)



**NBP1-72561C**

Goat anti-Mouse IgA/IgG/IgM (H+L) Secondary Antibody [DyLight 650] (Pre-adsorbed)

| Product Information         |                                                                                                                                                                                                                                                                                                                                            |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.5 ml                                                                                                                                                                                                                                                                                                                                     |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                   |
| Clonality                   | Polyclonal                                                                                                                                                                                                                                                                                                                                 |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                         |
| Isotype                     | IgG                                                                                                                                                                                                                                                                                                                                        |
| Conjugate                   | DyLight 650                                                                                                                                                                                                                                                                                                                                |
| Purity                      | Multi-step                                                                                                                                                                                                                                                                                                                                 |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                         |
| Product Description         |                                                                                                                                                                                                                                                                                                                                            |
| Host                        | Goat                                                                                                                                                                                                                                                                                                                                       |
| Species                     | Mouse                                                                                                                                                                                                                                                                                                                                      |
| Specificity/Sensitivity     | <p>This antibody was pre-adsorbed against Bovine, Horse, and Human Serum Proteins.</p> <p>No reaction was observed against bovine, horse or human serum proteins. This reagent is suitable for the detection of all mouse isotypes and chain combinations. ELISA was used to confirm specificity at less than 1% of the target signal.</p> |
| Immunogen                   | Mouse IgG IgA and IgM whole molecule                                                                                                                                                                                                                                                                                                       |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                                                                                                                          |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                            |
| Applications                | Western Blot, ELISA, Immunohistochemistry                                                                                                                                                                                                                                                                                                  |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry                                                                                                                                                                                                                                                                                                  |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                     |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-72561C](http://www.novusbio.com/reviews/submit/NBP1-72561C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

