

Product Datasheet

Goat anti-Mouse IgM Heavy Chain Secondary Antibody [DyLight 405] NBP1-72887

Unit Size: 0.1 mg

Store lyophilized antibody at 4C in the dark. Aliquot reconstituted liquid and store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-72887

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-72887



NBP1-72887**Goat anti-Mouse IgM Heavy Chain Secondary Antibody [DyLight 405]**

| Product Information | |
|------------------------------------|---|
| Unit Size | 0.1 mg |
| Concentration | LYOPH mg/ml |
| Storage | Store lyophilized antibody at 4C in the dark. Aliquot reconstituted liquid and store at -20C. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | 0.01% Sodium Azide |
| Reconstitution Instructions | Reconstitute with 100 ul deionized water (or equivalent). |
| Isotype | IgG |
| Conjugate | DyLight 405 |
| Purity | Multi-step |
| Buffer | Lyophilized from 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free |


| Product Description | |
|--------------------------------|--|
| Description | <p>Store vial at 4C prior to restoration. For extended storage aliquot contents and freeze at -20C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4C as an undiluted liquid. Dilute only prior to immediate use.</p> <p>This product was prepared from monospecific antiserum by immunoaffinity chromatography using Mouse IgM coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, Mouse IgM and Mouse Serum</p> |
| Host | Goat |
| Gene ID | 3507 |
| Gene Symbol | IGHM |
| Species | Mouse |
| Specificity/Sensitivity | No reaction was observed against other mouse heavy or light chain proteins. |
| Immunogen | Mouse IgM whole molecule |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details | |
|------------------------------------|---|
| Applications | Western Blot, Fluorophore-linked immunosorbent assay, Immunocytochemistry/Immunofluorescence |
| Recommended Dilutions | Western Blot 1:10000, Immunocytochemistry/ Immunofluorescence 1:5000, Fluorophore-linked immunosorbent assay 1:20000 |
| Application Notes | The emission spectra for this DyLight(TM) conjugate match the principle output wavelengths of most common fluorescence instrumentation. This product is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms. |



Images

Fluorophore-linked immunosorbent assay: Goat anti-Mouse IgM Heavy Chain Secondary Antibody [DyLight 405] [NBP1-72887] - Fluorescence absorption and emission spectra of DyLight 405.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-72887

| | |
|-------------|---------------------------------------|
| NB100-524 | NOD2 Antibody (2D9) - BSA Free |
| NBP1-96772 | Human IgM Heavy Chain Isotype Control |
| 7268-CT-100 | CTLA-4 [Unconjugated] |
| NB200-541 | Complement C4b/d Antibody (16D2) |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-72887

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

