

# Product Datasheet

## Goat anti-Guinea Pig IgG (H+L) Secondary Antibody [Alkaline Phosphatase] NBP1-74868

Unit Size: 1 mg

Store at 4C. After dilution, add an equal volume of glycerol to prevent loss of enzymatic activity and store at -20C. Avoid repeated freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-74868](http://www.novusbio.com/NBP1-74868)

Updated 10/30/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-74868](http://www.novusbio.com/reviews/destination/NBP1-74868)



**NBP1-74868**

Goat anti-Guinea Pig IgG (H+L) Secondary Antibody [Alkaline Phosphatase]

<b>Product Information</b>	
<b>Unit Size</b>	1 mg
<b>Concentration</b>	1.5 mg/ml
<b>Storage</b>	Store at 4C. After dilution, add an equal volume of glycerol to prevent loss of enzymatic activity and store at -20C. Avoid repeated freeze-thaw cycles.
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	0.05% Sodium Azide
<b>Isotype</b>	IgG
<b>Conjugate</b>	Alkaline Phosphatase
<b>Purity</b>	Affinity purified
<b>Buffer</b>	30 mM Triethanolamine, pH 7.2, 5 mM Magnesium Chloride, 0.1 mM Zinc Chloride, 1 % (w/v) BSA, Protease/IgG free
<b>Product Description</b>	
<b>Description</b>	Purity > 95% based on SDS-PAGE. Goat serum was obtained from healthy animals of US origin and under the care of a registered veterinarian.  Store undiluted liquid at 2-8C. For storage at -20C, dilute with an equal volume of glycerol to prevent loss of enzymatic activity. Prepare working dilutions prior to use and then discard.
<b>Host</b>	Goat
<b>Species</b>	Guinea Pig
<b>Reactivity Notes</b>	Based on IEP, no reactivity is observed to non-immunoglobulin guinea pig serum immunoglobulins
<b>Specificity/Sensitivity</b>	Based on IEP, this Goat anti-Guinea Pig IgG (H+L) Secondary Antibody [Alkaline Phosphatase] reacts with heavy gamma chains on guinea pig IgG and light chains on all guinea pig immunoglobulins.
<b>Immunogen</b>	This Goat anti-Guinea Pig IgG (H+L) Secondary Antibody [Alkaline Phosphatase] was developed against purified guinea pig IgG (H&L).
<b>Product Application Details</b>	
<b>Applications</b>	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin
<b>Recommended Dilutions</b>	Western Blot 1:200-1:5000, ELISA 1:200-1:5000, Immunohistochemistry 1:500-1:5000, Immunohistochemistry-Paraffin
<b>Application Notes</b>	This conjugate is suitable for all immunoassay applications.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-74868](http://www.novusbio.com/reviews/submit/NBP1-74868)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

