

Product Datasheet

Chicken anti-Mouse IgG (H+L) Secondary Antibody NBP1-75098

Unit Size: 2 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-75098

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-75098




NBP1-75098

Chicken anti-Mouse IgG (H+L) Secondary Antibody

Product Information	
Unit Size	2 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Product Description	
Host	Chicken
Species	Mouse
Reactivity Notes	Based on IEP, no reactivity is observed to non-immunoglobulin mouse serum immunoglobulins
Specificity/Sensitivity	Based on IEP, this Chicken anti-Mouse IgG (H+L) Secondary Antibody reacts with heavy gamma chains on mouse IgG and light chains on all mouse immunoglobulins.
Immunogen	This Chicken anti-Mouse IgG (H+L) Secondary Antibody was developed against purified mouse IgG (H&L).
Notes	Chicken products cannot be exported to Canada. Purity > 95% based on SDS-PAGE. Ig fraction was prepared using eggs obtained from healthy hens of US origin.
Product Application Details	
Applications	Western Blot, Chemiluminescence Immunoassay, ELISA, Flow Cytometry, Immunoassay, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunomicroscopy
Recommended Dilutions	Western Blot 1:100-1:2000, Flow Cytometry 1:10-1:1000, ELISA 1:100-1:2000, Immunohistochemistry 1:10 - 1:500, Immunocytochemistry/ Immunofluorescence 1:10 - 1:500, Immunoassay, Immunomicroscopy 1:10 - 1:500, Chemiluminescence Immunoassay
Application Notes	This antibody is suitable for all immunoassay applications.

Images

Immunocytochemistry/Immunofluorescence: Chicken anti-Mouse IgG (H +L) Secondary Antibody [NBP1-75098] - Analysis of Laminin in Vero cells. Cells were fixed in 4% formaldehyde for 20 minutes at room temperature. Image from verified customer review. Image using the DyLight 594 form of Chicken anti-Mouse IgG (H+L) Secondary Antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-75098

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

