

# Product Datasheet

## RHBDD3 Antibody [DyLight 594] NBP1-76250DL594

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-76250DL594](http://www.novusbio.com/NBP1-76250DL594)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-76250DL594](http://www.novusbio.com/reviews/destination/NBP1-76250DL594)



**NBP1-76250DL594**

RHBDD3 Antibody [DyLight 594]

| <b>Product Information</b>         |                                                                                                                                                                                                         |
|------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 ml                                                                                                                                                                                                  |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                             |
| <b>Storage</b>                     | Store at 4C in the dark.                                                                                                                                                                                |
| <b>Clonality</b>                   | Polyclonal                                                                                                                                                                                              |
| <b>Preservative</b>                | 0.05% Sodium Azide                                                                                                                                                                                      |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                     |
| <b>Conjugate</b>                   | DyLight 594                                                                                                                                                                                             |
| <b>Purity</b>                      | Peptide affinity purified                                                                                                                                                                               |
| <b>Buffer</b>                      | 50mM Sodium Borate                                                                                                                                                                                      |
| <b>Product Description</b>         |                                                                                                                                                                                                         |
| <b>Host</b>                        | Rabbit                                                                                                                                                                                                  |
| <b>Gene ID</b>                     | 25807                                                                                                                                                                                                   |
| <b>Gene Symbol</b>                 | RHBDD3                                                                                                                                                                                                  |
| <b>Species</b>                     | Human, Mouse, Rat                                                                                                                                                                                       |
| <b>Reactivity Notes</b>            | 0                                                                                                                                                                                                       |
| <b>Immunogen</b>                   | Antibody was raised against a 13 amino acid synthetic peptide from near the center of human RHBDD3. The immunogen is located within amino acids 170 - 220 of RHBDD3. Amino Acid Sequence: FRWLEPSERRLQV |
| <b>Notes</b>                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                       |
| <b>Product Application Details</b> |                                                                                                                                                                                                         |
| <b>Applications</b>                | Western Blot, ELISA, Immunocytochemistry/Immunofluorescence                                                                                                                                             |
| <b>Recommended Dilutions</b>       | Western Blot, ELISA, Immunocytochemistry/Immunofluorescence                                                                                                                                             |
| <b>Application Notes</b>           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                  |



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-76250DL594**

---

|                 |                                          |
|-----------------|------------------------------------------|
| NBP2-24891DL594 | Rabbit IgG Isotype Control [DyLight 594] |
| NBP1-92331PEP   | RHBDD3 Recombinant Protein Antigen       |
| NBL1-15342      | RHBDD3 Overexpression Lysate             |
| AF743           | VEGFR3/Flt-4 Antibody [Unconjugated]     |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-76250DL594](http://www.novusbio.com/reviews/submit/NBP1-76250DL594)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

