

# Product Datasheet

## ZNF346 Antibody [DyLight 755]

### NBP1-76540IR

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-76540IR](http://www.novusbio.com/NBP1-76540IR)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-76540IR](http://www.novusbio.com/reviews/destination/NBP1-76540IR)



**NBP1-76540IR**

ZNF346 Antibody [DyLight 755]

| Product Information         |                                                                                                                                                                                                                       |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                           |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                              |
| Clonality                   | Polyclonal                                                                                                                                                                                                            |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                    |
| Isotype                     | IgG                                                                                                                                                                                                                   |
| Conjugate                   | DyLight 755                                                                                                                                                                                                           |
| Purity                      | Peptide affinity purified                                                                                                                                                                                             |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                    |
| Product Description         |                                                                                                                                                                                                                       |
| Host                        | Rabbit                                                                                                                                                                                                                |
| Gene ID                     | 23567                                                                                                                                                                                                                 |
| Gene Symbol                 | ZNF346                                                                                                                                                                                                                |
| Species                     | Human, Mouse                                                                                                                                                                                                          |
| Reactivity Notes            | 0                                                                                                                                                                                                                     |
| Specificity/Sensitivity     | ZNF346 antibody is predicted to not cross-react with other ZNF family members. Three isoforms of ZNF346 are known to exist; this antibody will recognize the two longest isoforms                                     |
| Immunogen                   | Antibody was raised against a 20 amino acid synthetic peptide near the carboxy terminus of human ZNF346. The immunogen is located within the last 50 amino acids of ZNF346. Amino Acid Sequence: LGQIPMQRQPIQKDSTTLED |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                     |
| Product Application Details |                                                                                                                                                                                                                       |
| Applications                | Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                      |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin                                                                                                      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-76540IR**

---

|               |                                                  |
|---------------|--------------------------------------------------|
| NBP2-24891IR  | Rabbit IgG Isotype Control [DyLight 755]         |
| NBP2-22956    | Recombinant Human ZNF346 His Protein             |
| 203-IL-010    | IL-3 [Unconjugated]                              |
| H00023567-T01 | ZNF346 293T Cell Transient Overexpression Lysate |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-76540IR](http://www.novusbio.com/reviews/submit/NBP1-76540IR)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

