

# Product Datasheet

## Caspase-8 Antibody [DyLight 350] NBP1-76610UV

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-76610UV](http://www.novusbio.com/NBP1-76610UV)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-76610UV](http://www.novusbio.com/reviews/destination/NBP1-76610UV)



**NBP1-76610UV**

Caspase-8 Antibody [DyLight 350]

| Product Information         |                                                                                                                                                                                                                                        |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                 |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                            |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                               |
| Clonality                   | Polyclonal                                                                                                                                                                                                                             |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                     |
| Isotype                     | IgG                                                                                                                                                                                                                                    |
| Conjugate                   | DyLight 350                                                                                                                                                                                                                            |
| Purity                      | Peptide affinity purified                                                                                                                                                                                                              |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                     |
| Product Description         |                                                                                                                                                                                                                                        |
| Host                        | Rabbit                                                                                                                                                                                                                                 |
| Gene ID                     | 841                                                                                                                                                                                                                                    |
| Gene Symbol                 | CASP8                                                                                                                                                                                                                                  |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                                                      |
| Reactivity Notes            | 0                                                                                                                                                                                                                                      |
| Specificity/Sensitivity     | Depending on cell lines or tissues used, either full-length or other cleavage products may be observed.                                                                                                                                |
| Immunogen                   | Antibody was raised against a 16 amino acid synthetic peptide from near the carboxy terminus of human Caspase-8 isoform A. The immunogen is located within the last 50 amino acids of Caspase-8. Amino Acid Sequence: NYEVSNKDDKKNMGKQ |
| Notes                       | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                                                                                                                                      |
| Product Application Details |                                                                                                                                                                                                                                        |
| Applications                | Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                       |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin                                                                                                                       |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                 |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP1-76610UV**

---

|                      |                                             |
|----------------------|---------------------------------------------|
| NBP2-24891UV         | Rabbit IgG Isotype Control [DyLight 350]    |
| NBP2-58881PEP        | Caspase-8 Recombinant Protein Antigen       |
| 210-TA-005           | TNF-alpha [Unconjugated]                    |
| NBP2-54817-100Assays | Caspase-8 Activity Assay Kit (Colorimetric) |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-76610UV](http://www.novusbio.com/reviews/submit/NBP1-76610UV)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

