

Product Datasheet

IRAK3 Antibody - BSA Free

NBP1-76782

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-76782

Updated 10/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-76782



NBP1-76782

IRAK3 Antibody - BSA Free

Product Information

Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	70 kDa

Product Description


Description	Novus Biologicals Rabbit IRAK3 Antibody - BSA Free (NBP1-76782) is a polyclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	11213
Gene Symbol	IRAK3
Species	Human, Mouse, Rat
Specificity/Sensitivity	IRAK-M antibody is predicted to have no cross reactivity to other members of the IRAK family.
Immunogen	Antibody was raised against a peptide corresponding to 16 amino acids near the carboxy terminus of human IRAK-M. The immunogen is located within the last 50 amino acids of IRAK-M. Amino Acid Sequence: CSSKFSWDEYEYKKE

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1 ug/ml, ELISA 1:100-1:2000, Immunohistochemistry 2 ug/ml, Immunocytochemistry/ Immunofluorescence 10 ug/mL, Immunohistochemistry-Paraffin 2 ug/ml

Images

Western Blot: IRAK3 Antibody [NBP1-76782] - Analysis of IRAK3 in mouse spleen (M) and rat liver (R) tissue lysates with IRAK3 antibody using this antibody at 1 ug/ml.




Immunohistochemistry-Paraffin: IRAK3 Antibody [NBP1-76782] - Staining of human spleen tissue with antibody at 5 ug/mL.



Immunocytochemistry/ Immunofluorescence: IRAK3 Antibody - BSA Free [NBP1-76782] - Immunofluorescence of IRAK2 in A20 cells with IRAK2 antibody at 20 ug/mL.


Red: IRAK3 Antibody
Blue: DAPI staining



Immunohistochemistry: IRAK3 Antibody [NBP1-76782] - Rat Liver tissue with IRAK-M antibody at 10 ug/mL.



Immunohistochemistry: IRAK3 Antibody - BSA Free [NBP1-76782] - Immunohistochemical staining of rat liver tissue using IRAK3 antibody at 2 ug/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-76782

NBP1-76782PEP	IRAK3 Antibody Blocking Peptide
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-76782

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

