

Product Datasheet

SCF/c-kit Ligand Antibody [DyLight 594] NBP1-77030DL594

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-77030DL594

Updated 7/11/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-77030DL594



NBP1-77030DL594

SCF/c-kit Ligand Antibody [DyLight 594]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 594
Purity	Peptide affinity purified
Buffer	50mM Sodium Borate
Product Description	
Host	Rabbit
Gene ID	4254
Gene Symbol	KITLG
Species	Human, Mouse, Rat
Reactivity Notes	0
Specificity/Sensitivity	Human SCF/c-kit Ligand has 3 isoforms, including isoform 1 (273aa, 30.9kD), isoform 2 (245aa, 27.9kD) and isoform 3 (238aa, 26.7kD). Mouse SCF/c-kit Ligand has 2 isoforms, including isoform 1 (273aa, 30.6kD) and isoform 2 (245aa, 27.6kD). Rat SCF/c-kit Ligand also has 2 isoforms, including isoform 1 (273aa, 30.7kD) and isoform 2 (245aa, 27.7kD). NBP1-77030 can detect human and possibly detect mouse and rat.
Immunogen	Antibody was raised against a peptide corresponding to 18 amino acids near the center of human SCF/c-kit Ligand. The immunogen is located within amino acids 100-150 of SCF/c-kit Ligand. Amino Acid Sequence: KENSSKDLKKSFKSPEPR
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-77030DL594

NBP2-24891DL594	Rabbit IgG Isotype Control [DyLight 594]
NBP2-35256-100ug	Recombinant Human SCF/c-kit Ligand Protein
210-TA-005	TNF-alpha [Unconjugated]
255-SC-010	SCF/c-kit Ligand [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-77030DL594

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

