

# Product Datasheet

## SCF/c-kit Ligand Antibody [DyLight 350] NBP1-77030UV

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-77030UV](http://www.novusbio.com/NBP1-77030UV)

Updated 7/11/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-77030UV](http://www.novusbio.com/reviews/destination/NBP1-77030UV)



**NBP1-77030UV**

SCF/c-kit Ligand Antibody [DyLight 350]

<b>Product Information</b>	
<b>Unit Size</b>	0.1 ml
<b>Concentration</b>	Please see the vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Store at 4C in the dark.
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	0.05% Sodium Azide
<b>Isotype</b>	IgG
<b>Conjugate</b>	DyLight 350
<b>Purity</b>	Peptide affinity purified
<b>Buffer</b>	50mM Sodium Borate
<b>Product Description</b>	
<b>Host</b>	Rabbit
<b>Gene ID</b>	4254
<b>Gene Symbol</b>	KITLG
<b>Species</b>	Human, Mouse, Rat
<b>Reactivity Notes</b>	0
<b>Specificity/Sensitivity</b>	Human SCF/c-kit Ligand has 3 isoforms, including isoform 1 (273aa, 30.9kD), isoform 2 (245aa, 27.9kD) and isoform 3 (238aa, 26.7kD). Mouse SCF/c-kit Ligand has 2 isoforms, including isoform 1 (273aa, 30.6kD) and isoform 2 (245aa, 27.6kD). Rat SCF/c-kit Ligand also has 2 isoforms, including isoform 1 (273aa, 30.7kD) and isoform 2 (245aa, 27.7kD). NBP1-77030 can detect human and possibly detect mouse and rat.
<b>Immunogen</b>	Antibody was raised against a peptide corresponding to 18 amino acids near the center of human SCF/c-kit Ligand. The immunogen is located within amino acids 100-150 of SCF/c-kit Ligand. Amino Acid Sequence: KENSSKDLKKSFKSPEPR
<b>Notes</b>	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
<b>Product Application Details</b>	
<b>Applications</b>	Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
<b>Recommended Dilutions</b>	Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin
<b>Application Notes</b>	Optimal dilution of this antibody should be experimentally determined.



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-77030UV**

---

NBP2-24891UV	Rabbit IgG Isotype Control [DyLight 350]
NBP2-35256-100ug	Recombinant Human SCF/c-kit Ligand Protein
210-TA-005	TNF-alpha [Unconjugated]
255-SC-010	SCF/c-kit Ligand [Unconjugated]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-77030UV](http://www.novusbio.com/reviews/submit/NBP1-77030UV)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

