

# Product Datasheet

## Polypeptide GalNac Transferase 10/GALNT10 Antibody [PerCP] NBP1-77058PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-77058PCP](http://www.novusbio.com/NBP1-77058PCP)

Updated 10/23/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-77058PCP](http://www.novusbio.com/reviews/destination/NBP1-77058PCP)



**NBP1-77058PCP**

Polypeptide GalNAc Transferase 10/GALNT10 Antibody [PerCP]

| Product Information         |                                                                                                                                                                                                                      |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                               |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                          |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                             |
| Clonality                   | Polyclonal                                                                                                                                                                                                           |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                   |
| Isotype                     | IgG                                                                                                                                                                                                                  |
| Conjugate                   | PerCP                                                                                                                                                                                                                |
| Purity                      | Peptide affinity purified                                                                                                                                                                                            |
| Buffer                      | PBS                                                                                                                                                                                                                  |
| Product Description         |                                                                                                                                                                                                                      |
| Host                        | Rabbit                                                                                                                                                                                                               |
| Gene ID                     | 55568                                                                                                                                                                                                                |
| Gene Symbol                 | GALNT10                                                                                                                                                                                                              |
| Species                     | Human, Mouse, Rat                                                                                                                                                                                                    |
| Immunogen                   | Antibody was raised against an 18 amino acid synthetic peptide near the carboxy terminus of human GALNT10. The immunogen is located within amino acids 530 - 580 of GALNT10. Amino Acid Sequence: KGNQLWKYRKDKTLYHPV |
| Product Application Details |                                                                                                                                                                                                                      |
| Applications                | Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                     |
| Recommended Dilutions       | Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin                                                                                                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                               |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-77058PCP**

---

|               |                                                                       |
|---------------|-----------------------------------------------------------------------|
| NBP2-24891PCP | Rabbit IgG Isotype Control [PerCP]                                    |
| NBP1-85069PEP | Polypeptide GalNac Transferase 10/GALNT10 Recombinant Protein Antigen |
| 9077-GT-050   | GALNT14 [Unconjugated]                                                |
| 7575-GT-020   | Polypeptide GalNac Transferase 10/GALNT10                             |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-77058PCP](http://www.novusbio.com/reviews/submit/NBP1-77058PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

