

Product Datasheet

NLRP3/NALP3 Antibody - BSA Free NBP1-77080

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 2 Publications: 31

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-77080

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-77080




NBP1-77080

NLRP3/NALP3 Antibody - BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS
Target Molecular Weight	110 kDa
Product Description	
Host	Rabbit
Gene ID	114548
Gene Symbol	NLRP3
Species	Human, Mouse, E. coli
Reactivity Notes	E. coli reactivity reported in scientific literature (PMID: 31551961).
Immunogen	NLRP3/NALP3 Antibody was raised against a 16 amino acid synthetic peptide from near the amino terminus of human NALP3. The immunogen is located within amino acids 120 - 170 of NALP3.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation, SDS-Page
Recommended Dilutions	Western Blot 1 - 2 ug/mL, ELISA 1:100 - 1:2000, Immunohistochemistry 1:200, Immunocytochemistry/ Immunofluorescence 20 ug/mL, Immunoprecipitation reported in scientific literature (PMID 31551961), Immunohistochemistry-Paraffin 1:200, SDS-Page reported in scientific literature (PMID 31551961)

Images


Western Blot: NLRP3/NALP3 Antibody - BSA Free [NBP1-77080] - HMGB1 promoted the activation of canonical NLRP3 inflammasome in retinal IR injury. The intravitreal injection of HMGB1 inhibitor, GA, significantly suppressed mRNA (B) and protein (D) production of NLRP3. Data are shown as mean \pm SD. **P < 0.01 vs sham group. Image collected and cropped by CiteAb from the following publication (<https://www.jneuroinflammation.com/content/12/1/137>) licensed under a CC-BY license.




Immunocytochemistry/Immunofluorescence: NLRP3/NALP3 Antibody - BSA Free [NBP1-77080] - Staining of NALP3 in K562 cells with NALP3 antibody at 2 ug/mL.




Immunohistochemistry-Paraffin: NLRP3/NALP3 Antibody - BSA Free [NBP1-77080] - Analysis of NALP3 in mouse lung tissue using anti-NALP3 antibody. Image from verified customer review.



Western Blot: NLRP3/NALP3 Antibody - BSA Free [NBP1-77080] - Analysis of NALP3 in K562 cell lysate with NALP3 antibody at 1 ug/mL (A) in the absence and (B) in the presence of blocking peptide.




Immunocytochemistry/Immunofluorescence: NLRP3/NALP3 Antibody - BSA Free [NBP1-77080] - hADSC-Evs suppressed NLRP3 inflammasome activation in ocular surface of experimental mice under desiccating stress. Representative image of immunofluorescent staining of NLRP3 related genes in corneal and conjunctival epithelia of desiccating mice or treated with hADSC-EVs (*conjunctiva; open triangle cornea). Image collected and cropped by CiteAb from the following publication ([//pubmed.ncbi.nlm.nih.gov/32884023/](https://pubmed.ncbi.nlm.nih.gov/32884023/)) licensed under a CC-BY license.




hADSC-Evs suppressed NLRP3 inflammasome activation in ocular surface of experimental mice under desiccating stress. (A)

Representative image of immunofluorescent staining of NLRP3 related genes in corneal and conjunctival epithelia of desiccating mice or treated with hADSC-EVs (*conjunctiva; open triangle cornea). (B,C) qRT-PCR showed the mRNA levels of NLRP3 related genes in conjunctiva and cornea (four corneas/conjunctivas were mixed as a sample). Data was shown as mean \pm SEM. * $p < 0.05$.




Immunocytochemistry/ Immunofluorescence: NLRP3/NALP3 Antibody - BSA Free [NBP1-77080] - hADSC-Evs suppressed NLRP3

inflammasome activation in ocular surface of experimental mice under desiccating stress. (A) Representative image of immunofluorescent staining of NLRP3 related genes in corneal & conjunctival epithelia of desiccating mice or treated with hADSC-EVs (*conjunctiva; open triangle cornea). (B,C) qRT-PCR showed the mRNA levels of NLRP3 related genes in conjunctiva & cornea (four corneas/conjunctivas were mixed as a sample). Data was shown as mean \pm SEM. * $p < 0.05$. Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/32884023>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.




Western Blot: NLRP3/NALP3 Antibody - BSA Free [NBP1-77080] - HMGB1 promoted the activation of canonical NLRP3 inflammasome in retinal IR injury. The intravitreal injection of exogenous rHMGB1 significantly promoted the mRNA & protein production of NLRP3 (a, c) & ASC (e, g) & activation of caspase-1 (i). Thus, the intravitreal injection of HMGB1 inhibitor, GA, significantly suppressed the production of NLRP3 (b, d) & ASC (f, h) & activation of caspase-1 (j). Data are shown as mean \pm SD. * $P < 0.05$, ** $P < 0.01$ vs sham group Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/26224068>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.




Western Blot: NLRP3/NALP3 Antibody - BSA Free [NBP1-77080] - Co-immunoprecipitation assay: Figure showing interaction of NLRP3, NLRC4 & NAIP with ASC during CFT073 infection in THP-1m.

Cytoplasmic extracts of THP-1m infected with or without CFT073, 6 h post stimulation were subjected to IP & WB by various antibodies (ASC, NLRP3, NLRC4, & NAIP) as described in material & methods. Cytoplasmic extract used as input was subjected to check NLRP3, ASC, NLRC4 & NAIP expression as shown in figure. ASC antibody was used to co-immunoprecipitate the complex from cytoplasmic extracts & then immunoblotted to detect NLRP3, NLRC4 & NAIP. Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/31551961>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Western Blot: NLRP3/NALP3 Antibody - BSA Free [NBP1-77080] - HMGB1 regulated the activation of NLRP3, rather than the activation of caspase-8, via NF- κ B pathway. a, b NF- κ B was activated in ischemic retina at the early stage after reperfusion. c, d Intraocular injection rHMGB1 or GA could promote or suppress the production of phosphor-NF- κ B p65. e, f Intravitreal injection of NF- κ B p65 inhibitor, JSH-23, significantly reduced the activation of NLRP3 & decreased the processing of IL-1 β , g rather than caspase-8. However, intravitreal injection of caspase-8 inhibitor, Z-IETD-fmk, decreased the production of phosphor-NF- κ B p65, obviously (h). Data are shown as mean \pm SD. *P < 0.05, **P < 0.01 vs sham group Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/26224068>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

Chi W, Hua X, Chen X et al. Mitochondrial DNA oxidation induces imbalanced activity of NLRP3/NLRP6 inflammasomes by activation of caspase-8 and BRCC36 in dry eye. *J. Autoimmun.* 2017-02-23 [PMID: 28238526]

Schultz V, Barrie JA, Donald CL et al. Oligodendrocytes are susceptible to Zika virus infection in a mouse model of perinatal exposure: Implications for CNS complications *Glia* 2021-08-01 [PMID: 33942402]

Chen XC, Wu D, Wu HL et al. Metformin improves renal injury of MRL/lpr lupus-prone mice via the AMPK/STAT3 pathway *Lupus Science & Medicine* 2022-04-11 [PMID: 35414608] (In vivo assay)

Wu F, Wei H, Hu Y et al. Upregulation of P2X7 Exacerbates Myocardial Ischemia-Reperfusion Injury through Enhancing Inflammation and Apoptosis in Diabetic Mice *Journal of immunology (Baltimore, Md. : 1950)* 2023-05-05 [PMID: 37144844] (WB, Mouse)

Details:

1:1000 WB dilution

Yu Z, Yazdanpanah G, Alam J et al. Induction of Innate Inflammatory Pathways in the Corneal Epithelium in the Desiccating Stress Dry Eye Model *Investigative ophthalmology & visual science* 2023-04-03 [PMID: 37036417] (WB, Mouse)

Gupta A, Singh K, Fatima S et al. Neutrophil Extracellular Traps Promote NLRP3 Inflammasome Activation and Glomerular Endothelial Dysfunction in Diabetic Kidney Disease *Nutrients* 2022-07-20 [PMID: 35889923] (WB, Mouse)

Details:

1:1000 WB dilution

Jin L, Jin F, Guo S et al. Metformin Inhibits NLR Family Pyrin Domain Containing 3 (NLRP)-Relevant Neuroinflammation via an Adenosine-5-Monophosphate-Activated Protein Kinase (AMPK)-Dependent Pathway to Alleviate Early Brain Injury After Subarachnoid Hemorrhage in Mice *Frontiers in Pharmacology* 2022-03-17 [PMID: 35370693] (IF/IHC, Mouse)

Yu, C, Chen, P Et al. hADSCs derived extracellular vesicles inhibit NLRP3 inflammasome activation and dry eye. *Sci Rep* 2020-09-03 [PMID: 32884023] (WB, Mouse)

Gaul S, Schaeffer Km, Opitz L Et Al. Extracellular NLRP3 inflammasome particles are internalized by human coronary artery smooth muscle cells and induce pro-atherogenic effects *Scientific reports* 2021-07-26 [PMID: 34312415] (WB)

Hua X, Yuan X et al. Desiccating stress worsens alkali burn injury by magnifying caspase-8-induced imbalance of NLRP3 and NLRP6. *J Allergy Clin Immunol* 2017-01-10 [PMID: 28479336] (IF/IHC, Mouse)

Verma V, Kumar P, Gupta S et al. Alpha-Hemolysin of uropathogenic *E. coli* regulates NLRP3 inflammasome activation and mitochondrial dysfunction in THP-1 macrophages *Sci Rep* 2020-07-28 [PMID: 32724079] (IP, WB, Human)

Toyohara T, Roudnicky F, Florido MHC et al. Patient hiPSCs Identify Vascular Smooth Muscle Arylacetamide Deacetylase as Protective against Atherosclerosis *Cell Stem Cell* 2020-05-14 [PMID: 32413331] (IHC-P, Mouse)

More publications at <http://www.novusbio.com/NBP1-77080>





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-77080

NBP1-77080PEP	NLRP3/NALP3 Antibody Blocking Peptide
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-77080

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

