

Product Datasheet

ZEB1 Antibody - BSA Free NBP1-77178

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-77178

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-77178



NBP1-77178

ZEB1 Antibody - BSA Free

Product Information

Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	130 kDa

Product Description


Host	Rabbit
Gene ID	6935
Gene Symbol	ZEB1
Species	Human
Marker	Mesenchymal Cells Marker
Specificity/Sensitivity	Multiple isoforms of ZEB1 are known to exist.
Immunogen	Antibody was raised against a 15 amino acid synthetic peptide near the center of human ZEB1. The immunogen is located within amino acids 530 - 580 of ZEB1. Amino Acid Sequence: HYDLKQPTQPPPLPA

Product Application Details


Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1 ug/ml, ELISA 1:100-1:2000, Immunohistochemistry 5 ug/ml, Immunocytochemistry/ Immunofluorescence 20 ug/mL, Immunohistochemistry-Paraffin 5 ug/ml

Images


Western Blot: ZEB1 Antibody [NBP1-77178] - HeLa cell lysate with ZEB1 antibody at 1ug/mL.




Immunocytochemistry/Immunofluorescence: ZEB1 Antibody [NBP1-77178] - Staining of ZEB1 in K562 cells with ZEB1 antibody at 2.5 ug/ml.




Immunohistochemistry-Paraffin: ZEB1 Antibody [NBP1-77178] - Staining of human liver tissue with antibody at 5 ug/mL.



Immunocytochemistry/Immunofluorescence: ZEB1 Antibody [NBP1-77178] - Staining of ZEB1 in HeLa cells with ZEB1 antibody at 20 ug/mL.



Immunocytochemistry/Immunofluorescence: ZEB1 Antibody [NBP1-77178] - Hela cells at 20 ug/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-77178

NBP1-77178PEP	ZEB1 Antibody Blocking Peptide
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-77178

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

