

Product Datasheet

cIAP-1/HiAP-2 Antibody [DyLight 755] NBP1-77193IR

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-77193IR

Updated 7/11/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-77193IR



NBP1-77193IR

cIAP-1/HiAP-2 Antibody [DyLight 755]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 755
Purity	Peptide affinity purified
Buffer	50mM Sodium Borate
Product Description	
Host	Rabbit
Gene ID	329
Gene Symbol	BIRC2
Species	Human, Mouse, Rat
Reactivity Notes	0
Specificity/Sensitivity	Human cIAP-1/HiAP-2 has 2 isoforms, including isoform 1 (618aa, 70kD) and isoform 2 (569aa, 64kD). Mouse cIAP-1/HiAP-2 has one isoform (612aa, 70kD) and Rat cIAP-1/HiAP-2 also has one isoform (589aa, 67kD). NBP1-77193 can detect human, mouse and rat.
Immunogen	Antibody was raised against a peptide corresponding to 14 amino acids near the carboxy terminus of human cIAP-1/HiAP-2. c-IAP antibody detects both c-IAP1 and c-IAP2. The immunogen is located within amino acids 540-590 of cIAP-1/HiAP-2. Amino Acid Sequence: SLEEQLRRLQEERT
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin
Application Notes	Optimal dilution of this antibody should be experimentally determined.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-77193IR

NBP2-24891IR	Rabbit IgG Isotype Control [DyLight 755]
NBP1-90133PEP	cIAP-1/HiAP-2 Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]
818-IA-050	cIAP-1/HiAP-2

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-77193IR

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

