

Product Datasheet

R-Spondin 1 Antibody - BSA Free NBP1-77354

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-77354

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-77354



NBP1-77354

R-Spondin 1 Antibody - BSA Free


Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	29 kDa

Product Description	
Host	Rabbit
Gene ID	284654
Gene Symbol	RSPO1
Species	Human, Mouse
Specificity/Sensitivity	At least three isoforms of R-Spondin 1 are known to exist; this antibody will detect the largest and smallest isoforms. R-Spondin 1 antibody will not cross-react with other RSPO family members. Human R-Spondin 1 has 3 isoforms, including isoform 1 (263aa, 29kD), isoform 2 (236aa, 26kD), and isoform 3 (200aa, 22kD). Mouse R-Spondin 1 had 1 isoform (265aa, 29kD) and rat R-Spondin 1 had 1 isoform (262aa, 29kD). NBP1-77354 can detect human, mouse and rat.
Immunogen	Antibody was raised against a peptide corresponding to 16 amino acids near the amino terminus of human R-Spondin 1. The immunogen is located within the first 50 amino acids of R-Spondin 1. Amino Acid Sequence: RGIKGRQRRISAEGS


Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 0.5-2 ug/mL, ELISA 1:100-1:2000, Immunohistochemistry 2.5-5 ug/mL, Immunocytochemistry/ Immunofluorescence 5 ug/ml, Immunohistochemistry-Paraffin 5 ug/ml

Images


Western Blot: R-Spondin 1 Antibody [NBP1-77354] - Loading: 15 ug of spleen lysates per lane. Antibodies: RSPO1, 5171 (Lane 1: 0.5 ug/mL, Lane 2: 1 ug/mL), 1h incubation at RT in 5% NFDN/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.




Immunocytochemistry/Immunofluorescence: R-Spondin 1 Antibody - BSA Free [NBP1-77354] - R-Spondin 1 Antibody [NBP1-77354] - Staining of RSPO1 in human spleen tissue with RSPO1 antibody at 5 ug/ml.




Immunohistochemistry: R-Spondin 1 Antibody [NBP1-77354] - RSPO1 in human spleen tissue with RSPO1 antibody at 5 ug/ml.




Western Blot: R-Spondin 1 Antibody [NBP1-77354] - Analysis in human heart tissue lysate with antibody at (A) 1 and (B) 2 ug/mL.




Western Blot: R-Spondin 1 Antibody [NBP1-77354] - Loading: 30 ng of lysates per lane. Lane 1: 0.5ug/ml, Lane 2: 1ug/ml




Western Blot: R-Spondin 1 Antibody [NBP1-77354] - Loading: 30 ng of human RSPO1 recombinant protein per lane. Antibodies: RSPO1, 5171 (Lane 1: 0.5 ug/mL, Lane 2: 1 ug/mL, Lane 3: 2 ug/mL), 1h incubation at RT in 5% NFDN/TBST. Secondary: Goat anti-rabbit IgG HRP conjugate at 1:10000 dilution.



Immunohistochemistry-Paraffin: R-Spondin 1 Antibody [NBP1-77354] - Human spleen tissue with RSPO1 antibody at 2.5 ug/ml.



Immunohistochemistry: R-Spondin 1 Antibody [NBP1-77354] - Analysis of RSPO1 in human heart tissue with RSPO1 antibody at 2.5 ug/ml.



Publications

MY Turco, L Gardner, RG Kay, RS Hamilton, M Prater, MS Hollinshea, A McWhinnie, L Esposito, R Fernando, H Skelton, F Reimann, FM Gribble, A Sharkey, SGE Marsh, S O'Rahilly, M Hemberger, GJ Burton, A Moffett
Trophoblast organoids as a model for maternal-fetal interactions during human placentation *Nature*, 2018-11-28;0(0):. 2018-11-28 [PMID: 30487605]

Takahashi Y, Sato S, Kurashima Y et al. A Refined Culture System for Human Induced Pluripotent Stem Cell-Derived Intestinal Epithelial Organoids. *Stem Cell Reports* 2017-12-01 [PMID: 29233552] (Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-77354

NBP1-77354PEP	R-Spondin 1 Antibody Blocking Peptide
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-77354

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

