

Product Datasheet

R-Spondin 1 Antibody [DyLight 755] NBP1-77354IR

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-77354IR

Updated 7/11/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-77354IR



NBP1-77354IR

R-Spondin 1 Antibody [DyLight 755]

| Product Information | |
|------------------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Polyclonal |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG |
| Conjugate | DyLight 755 |
| Purity | Peptide affinity purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Host | Rabbit |
| Gene ID | 284654 |
| Gene Symbol | RSPO1 |
| Species | Human, Mouse |
| Reactivity Notes | 0 |
| Specificity/Sensitivity | At least three isoforms of R-Spondin 1 are known to exist; this antibody will detect the largest and smallest isoforms. R-Spondin 1 antibody will not cross-react with other RSPO family members. Human R-Spondin 1 has 3 isoforms, including isoform 1 (263aa, 29kD), isoform 2 (236aa, 26kD), and isoform 3 (200aa, 22kD). Mouse R-Spondin 1 had 1 isoform (265aa, 29kD) and rat R-Spondin 1 had 1 isoform (262aa, 29kD). NBP1-77354 can detect human, mouse and rat. |
| Immunogen | Antibody was raised against a peptide corresponding to 16 amino acids near the amino terminus of human R-Spondin 1. The immunogen is located within the first 50 amino acids of R-Spondin 1. Amino Acid Sequence: RGIK GKRRISAEGS |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |
| Product Application Details | |
| Applications | Western Blot, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin |
| Recommended Dilutions | Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-77354IR

| | |
|----------------|--|
| NBP2-24891IR | Rabbit IgG Isotype Control [DyLight 755] |
| NBP3-21297PEP | R-Spondin 1 Recombinant Protein Antigen |
| 5036-WN-010 | Wnt-3a [Unconjugated] |
| 7150-RS-025/CF | R-Spondin 1 [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-77354IR

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

