

Product Datasheet

WDR37 Antibody - BSA Free

NBP1-81449

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-81449

Updated 4/22/2026 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-81449



NBP1-81449

WDR37 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Host	Rabbit
Gene ID	22884
Gene Symbol	WDR37
Species	Human
Reactivity Notes	Rat (81%).
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: TARQTKQKRKSHLSIRRTNSSEQERTGLPRDMLEGQDSKLPSSVRSTLLELF GQIEREFENLYIENLELRREIDTLNERLAAEG


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:20 - 1:50, Immunohistochemistry-Paraffin 1:20 - 1:50
Application Notes	IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Immunohistochemistry-Paraffin: WDR37 Antibody [NBP1-81449] - Staining of human lymph node shows strong cytoplasmic positivity in non-germinal center cells.




Analysis in human cell lines RT-4 and U2OS using Anti-WDR37 antibody. Corresponding WDR37 RNA-seq data are presented for the same cell lines. Loading control: Anti-PFN1.




Immunohistochemistry-Paraffin: WDR37 Antibody [NBP1-81449] - Staining of human cerebral cortex shows strong positivity in neuropil.




Immunohistochemistry-Paraffin: WDR37 Antibody [NBP1-81449] - Staining of human cerebral cortex, kidney, liver and lymph node using Anti-WDR37 antibody NBP1-81449 (A) shows similar protein distribution across tissues to independent antibody NBP1-81450 (B).



Immunohistochemistry-Paraffin: WDR37 Antibody [NBP1-81449] - Staining of human kidney shows weak cytoplasmic positivity in cells in tubules.



Immunohistochemistry-Paraffin: WDR37 Antibody [NBP1-81449] - Staining of human liver shows very weak positivity in hepatocytes as expected.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-81449

NBP1-81449PEP	WDR37 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-81449

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

