

Product Datasheet

TMEM174 Antibody - BSA Free NBP1-82023

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-82023

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-82023



NBP1-82023

TMEM174 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol


Product Description	
Description	Novus Biologicals Rabbit TMEM174 Antibody - BSA Free (NBP1-82023) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-TMEM174 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	134288
Gene Symbol	TMEM174
Species	Human
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: LITSGGAAAAMSSPPQYYTIYPQDNSAFVVDEGCLSFTDGGNHRPNPDVDQLE ETQLEEEACACFSPPPYEEIYSLPR

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry 1:1000 - 1:2500, Immunohistochemistry-Paraffin 1:1000 - 1:2500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.




Images


Immunohistochemistry-Paraffin: TMEM174 Antibody [NBP1-82023] - Staining of human lymph node shows low expression as expected.



Immunohistochemistry-Paraffin: TMEM174 Antibody [NBP1-82023] - Staining of human kidney shows high expression.



Immunohistochemistry-Paraffin: TMEM174 Antibody [NBP1-82023] - Staining in human kidney and lymph node tissues using anti-TMEM174 antibody. Corresponding TMEM174 RNA-seq data are presented for the same tissues.



Publications

Miyazaki-Anzai S, Keenan A, Blaine J, Miyazaki M Targeted Disruption of a Proximal Tubule Specific TMEM174 Gene in Mice Causes Hyperphosphatemia and Vascular Calcification Journal of the American Society of Nephrology : JASN 2022-04-22 [PMID: 35459732] (ICC/IF, IF/IHC, WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-82023

NBP1-82023PEP	TMEM174 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-82023

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

