

Product Datasheet

HADHA Antibody - BSA Free

NBP1-83387

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-83387

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-83387



NBP1-83387

HADHA Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Host	Rabbit
Gene ID	3030
Gene Symbol	HADHA
Species	Human, Mouse
Reactivity Notes	Rat (85%).
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: IEYLEEVAITFAKGLADKKISPKRDKGLVEKLTAYAMTIPFVRQQVYKKVEEKVR KQTKGLYPAPLKIIDVVKTGIEQGS DAGYLCESQKFGELVMTKESKALMGLYHG QVLCKKNKFGAPQKDVKHLAILGAGLMGAGIAQVSVDK

Product Application Details


Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Knockdown Validated
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:200 - 1:500, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:200 - 1:500, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.

Images


Western Blot: HADHA Antibody [NBP1-83387] - Analysis in Hep-G2 cells transfected with control siRNA, target specific siRNA probe #1, using Anti-HADHA antibody. Remaining relative intensity is presented. Loading control: Anti-GAPDH.




Immunocytochemistry/Immunofluorescence: HADHA Antibody [NBP1-83387] - Staining of human cell line U-2 OS shows localization to mitochondria. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: HADHA Antibody [NBP1-83387] - Staining of human small intestine shows strong granular cytoplasmic positivity in glandular cells.




Western Blot: HADHA Antibody [NBP1-83387] - Analysis in mouse cell line NIH-3T3 and rat cell line NBT-II.




Immunohistochemistry-Paraffin: HADHA Antibody [NBP1-83387] - Staining of human kidney shows moderate to strong cytoplasmic positivity in cells in tubules.



Immunohistochemistry-Paraffin: HADHA Antibody [NBP1-83387] - Staining of human liver shows strong cytoplasmic positivity in hepatocytes.



Immunohistochemistry-Paraffin: HADHA Antibody [NBP1-83387] - Staining of human pancreas shows moderate cytoplasmic positivity in exocrine glandular cells.



Publications

Maria Weiss, Sara Hettrich, Theresa Hofmann, Salma Hachim, Stefan Günther, Thomas Braun, Thomas Boettger
Mitolnc controls cardiac BCAA metabolism and heart hypertrophy by allosteric activation of BCKDH Nucleic Acids
Research 2024-06-24 [PMID: 38567728]

Gerin I, Clerbaux LA, Haumont O et al. Expression of miR-33 from an SREBP2 intron inhibits cholesterol export and
fatty acid oxidation. J Biol Chem 2010-10-01 [PMID: 20732877]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-83387

NBP1-83387PEP	HADHA Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-83387

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

