

Product Datasheet

DGCR2 Recombinant Protein Antigen NBP1-84255PEP

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-84255PEP

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-84255PEP



NBP1-84255PEP**DGCR2 Recombinant Protein Antigen**

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Preservative	No Preservative
Purity	>80% by SDS-PAGE and Coomassie blue staining
Buffer	PBS and 1M Urea, pH 7.4.
Target Molecular Weight	39 kDa

Product Description	
Description	A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human DGCR2. Source: <i>E. coli</i> Amino Acid Sequence: RFSRKCP TGWHHYEGTASCYRVYLSGENYWDAAQTCQRLNGSLATFSTDQE LRFVLAQEWDQPER SFGWKDQRKLWVG YQYVITGRNRSLEGRWEVAFKGSS EVFLPPDPIFASAMSENDNVFCAQLQCFHFPTLRHHD LHSWHAESCYEKSSFL CKRSQTCVDIKDNVVDEGFYFTPK Fusion Tag: N-terminal His6ABP (ABP = Albumin Binding Protein derived from Streptococcal Protein G) This product is intended to be used as a blocking antigen for antibody competition assays. Any other use of this antigen is done at the risk of the user. The use of this product for commercial production is strictly prohibited. Please contact technical support if you have any questions.
Gene ID	9993
Gene Symbol	DGCR2
Species	Human

Product Application Details	
Applications	Antibody Competition
Recommended Dilutions	Antibody Competition 10 - 100 molar excess



Application Notes	This recombinant antigen is only intended to be used as a blocking agent to confirm antibody specificity with the corresponding antibody, catalog number NBP1-84255. It is purified by IMAC chromatography, and the expected concentration is greater than 0.5 mg/ml. For current lot information, including availability, please contact our technical support team click nb-technical@bio-techne.com For further blocking peptide related information and a protocol, click here.
--------------------------	--





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-84255PEP

NBP2-56538PEP	DGCR2 Recombinant Protein Antigen
202-IL-010	IL-2 [Unconjugated]
MAB10443-100	DGCR2 Antibody (1019908) [Unconjugated]
291-G1-200	IGF-I/IGF-1 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-84255PEP

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

