

Product Datasheet

5'-Nucleotidase/CD73 Antibody - BSA Free NBP1-85740

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 4 Publications: 12

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-85740

Updated 2/21/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-85740



NBP1-85740

5'-Nucleotidase/CD73 Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Host	Rabbit
Gene ID	4907
Gene Symbol	NT5E
Species	Human, Porcine
Reactivity Notes	Porcine reactivity reported from a verified customer review.
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: EFDNGVEGLIEPLLKEAKFPILSANIKAKGPLASQISGLYLPYKVLPVGDEVVGIV GYTSKETPFLSNPGTNLVFEDEITALQPEVDKCLKTLNVNKKIIALGHSGFEMDKLIA QK


Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, ELISA Reported in scientific literature (PMID:33488292)., Immunohistochemistry 1:500 - 1:1000, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:500 - 1:1000
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization, Use PFA/Triton X-100.

Images


Western Blot: 5'-Nucleotidase/CD73 Antibody [NBP1-85740] - Analysis in human cell line U-87 MG.




Immunocytochemistry/Immunofluorescence: 5'-Nucleotidase/CD73 Antibody [NBP1-85740] - Staining for human aortic valve endothelial cells. Image from a verified customer review.




Immunohistochemistry-Paraffin: 5'-Nucleotidase/CD73 Antibody [NBP1-85740] - Staining of human liver shows strong membranous positivity in hepatocytes.




Immunocytochemistry/Immunofluorescence: 5'-Nucleotidase/CD73 Antibody [NBP1-85740] - Staining of human cell line LHCN-M2 shows localization to nucleoplasm, plasma membrane & cytosol. Antibody staining is shown in green.




Immunohistochemistry: 5'-Nucleotidase/CD73 Antibody [NBP1-85740] - Staining of human aortic valve staining for CD73. Image from verified customer review.




Immunohistochemistry: 5'-Nucleotidase/CD73 Antibody [NBP1-85740] - Staining for pig aortic valve. Image from verified customer review.




Immunohistochemistry-Paraffin: 5'-Nucleotidase/CD73 Antibody [NBP1-85740] - Staining of human endometrium shows strong membranous positivity in glandular cells.



Immunohistochemistry-Paraffin: 5'-Nucleotidase/CD73 Antibody [NBP1-85740] - Staining of human duodenum shows strong positivity in apical membrane in glandular cells.



Immunohistochemistry-Paraffin: 5'-Nucleotidase/CD73 Antibody [NBP1-85740] - Staining of human endometrium shows strong positivity in apical membrane in glandular cells.



Publications

Ahmed MS, Farag AB, Boys IN et al. FDA approved drugs with antiviral activity against SARS-CoV-2: From structure-based repurposing to host-specific mechanisms *Biomedicine & Pharmacotherapy* 2023-06-01 [PMID: 37068330] (Western Blot)

Monguio-Tortajada M, Prat-Vidal C, Martinez-Falguera D et al. Acellular cardiac scaffolds enriched with MSC-derived extracellular vesicles limit ventricular remodelling and exert local and systemic immunomodulation in a myocardial infarction porcine model *Theranostics* 2022-06-06 [PMID: 35832072] (FLOW, IHC-Fr, Porcine)

Pang L, Ng KT, Liu J et al. Plasmacytoid dendritic cells recruited by HIF-1 alpha /eADO/ADORA1 signaling induce immunosuppression in hepatocellular carcinoma *Cancer letters* 2021-09-16 [PMID: 34536555] (IF/IHC, Human)

MuNOz-Godlnez R, de Lourdes Mora-Garcia M, Weiss-Steider B et al. Detection of CD39 and a Highly Glycosylated Isoform of Soluble CD73 in the Plasma of Patients with Cervical Cancer: Correlation with Disease Progression Mediators of Inflammation 2020-12-07 [PMID: 33488292] (ELISA, WB, Human)

Details:

Platelet-free plasma from human blood underwent ELISA.

Gong W, Zhang S, Zong Y et al. HPV-16 Infection Is Associated with a High Content of CD39 and CD73 Ectonucleotidases in Cervical Samples from Patients with CIN-1 Mediators of Inflammation 2019-05-07 [PMID: 31205451] (IF/IHC, Human)

Calabrese G, Giuffrida R, Forte S et al. Human adipose-derived mesenchymal stem cells seeded into a collagen-hydroxyapatite scaffold promote bone augmentation after implantation in the mouse *Sci Rep* 2017-08-02 [PMID: 28769083] (ICC/IF, Human)

Bullen JW, Tchernyshyov I, Holewinski RJ et al. Protein kinase A-dependent phosphorylation stimulates the transcriptional activity of hypoxia-inducible factor 1 *Sci Signal* 2016-05-31 [PMID: 27245613] (WB, Human)

Thibaudin M, Chaix M, Boidot R et al. Human ectonucleotidase-expressing CD25high Th17 cells accumulate in breast cancer tumors and exert immunosuppressive functions. *Oncoimmunology* 2015-07-06 [PMID: 26942062] (WB)

Serra S, Horenstein AL, Vaisitti T et al. CD73-generated extracellular adenosine in chronic lymphocytic leukemia creates local conditions counteracting drug-induced cell death. *Blood* 2011-12-01 [PMID: 21998208]

D'Alimonte I, Nargi E, Zuccarini M et al. Potentiation of temozolomide antitumor effect by purine receptor ligands able to restrain the in vitro growth of human glioblastoma stem cells. *Purinergic Signal*. 2015-05-15 [PMID: 25976165] (WB, Human)

Alberston P, Dex S, Popov C et al. Loss of tenomodulin results in reduced self-renewal and augmented senescence of tendon stem/progenitor cells. *Stem Cells Dev*. 2014-10-28 [PMID: 25351164]

Snider Nt, Altshuler Pj, Wan S et al. Alternative splicing of human nT5E in cirrhosis and hepatocellular carcinoma produces a negative regulator of ecto-5'-nucleotidase (CD73). *Mol. Biol. Cell*. 2014-10-08 [PMID: 25298403]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-85740

NBP1-85740PEP	5'-Nucleotidase/CD73 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-85740

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

