

# Product Datasheet

## ZNF313 Antibody - BSA Free

### NBP1-87204

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

#### Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-87204](http://www.novusbio.com/NBP1-87204)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-87204](http://www.novusbio.com/reviews/destination/NBP1-87204)



**NBP1-87204**

ZNF313 Antibody - BSA Free

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                           |
| Clonality           | Polyclonal                                                                                                       |
| Preservative        | 0.02% Sodium Azide                                                                                               |
| Isotype             | IgG                                                                                                              |
| Purity              | Immunogen affinity purified                                                                                      |
| Buffer              | PBS (pH 7.2) and 40% Glycerol                                                                                    |


| Product Description |                                                                                                                                                                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | Novus Biologicals Rabbit ZNF313 Antibody - BSA Free (NBP1-87204) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-ZNF313 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host                | Rabbit                                                                                                                                                                                                                                               |
| Gene ID             | 55905                                                                                                                                                                                                                                                |
| Gene Symbol         | RNF114                                                                                                                                                                                                                                               |
| Species             | Human, Mouse, Rat                                                                                                                                                                                                                                    |
| Reactivity Notes    | Reactivity reported in scientific literature (PMID: 23435261)                                                                                                                                                                                        |
| Immunogen           | This antibody was developed against Recombinant Protein corresponding to amino acids:<br>SALAPGVRAVELERQIESTETSCHGCRKNFFLSKIRSHVATCSKYQNYIMEGVK<br>ATIKDASLQPRNVNRYTFPCPYCPEKNFDQEGLVEHCKLFHSTDTKSVVCPIC<br>ASMPW                                    |

| Product Application Details |                                                                                                                                                                    |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry                                                         |
| Recommended Dilutions       | Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:200 - 1:500, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:200 - 1:500 |
| Application Notes           | For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.                                                      |




## Images

Western Blot: ZNF313 Antibody [NBP1-87204] - Lane 1: NIH-3T3 cell lysate (Mouse embryonic fibroblast cells). Lane 2: NBT-II cell lysate (Rat Wistar bladder tumor cells).




Immunocytochemistry/Immunofluorescence: ZNF313 Antibody [NBP1-87204] - Staining of human cell line A-431 shows localization to plasma membrane & cytosol. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: ZNF313 Antibody [NBP1-87204] - Staining of human kidney shows strong cytoplasmic positivity in cells in tubules.




Western Blot: ZNF313 Antibody [NBP1-87204] - Lane 1: Marker [kDa] 230, 130, 95, 72, 56, 36, 28, 17, 11. Lane 2: Human cell line RT-4. Lane 3: Human cell line U-251MG sp



Immunohistochemistry-Paraffin: ZNF313 Antibody [NBP1-87204] - Staining of human kidney shows strong cytoplasmic positivity in cells in tubules.




Immunohistochemistry-Paraffin: ZNF313 Antibody [NBP1-87204] - Staining of human testis shows strong cytoplasmic and nuclear positivity in cells in seminiferous ducts.



Immunohistochemistry-Paraffin: ZNF313 Antibody [NBP1-87204] - Staining of human duodenum shows moderate cytoplasmic positivity in glandular cells.



Immunohistochemistry-Paraffin: ZNF313 Antibody [NBP1-87204] - Staining of human liver shows moderate cytoplasmic positivity in hepatocytes.



## Publications

Stadler C, Rexhepaj E, Singan VR et al. Immunofluorescence and fluorescent-protein tagging show high correlation for protein localization in mammalian cells. Nat Methods 2013-04-01 [PMID: 23435261]



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP1-87204**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| NBP1-87204PEP | ZNF313 Recombinant Protein Antigen                  |
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-87204](http://www.novusbio.com/reviews/submit/NBP1-87204)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

