

Product Datasheet

PUM3 Antibody - BSA Free

NBP1-87296

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-87296

Updated 12/2/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-87296



NBP1-87296

PUM3 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol


Product Description	
Description	Novus Biologicals Rabbit PUM3 Antibody - BSA Free (NBP1-87296) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	9933
Gene Symbol	KIAA0020
Species	Human
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: Mouse (85%), Rat (85%)
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: EHAQEVVLDKSACVLVSDILGSATGDVQPTMNAIASLAATGLHPGGKDGEHIA EHPAGHLVLKWLIEQDKMKENGREGCFAKTLVEHVGMKNLKSASVNRGAI LSSLLQSCDLEVANKVKAALKSLIPTLEKTKSTSKGI

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:20 - 1:50, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:20 - 1:50
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. Immunocytochemistry/Immunofluorescence Fixation Permeabilization: Use PFA/Triton X-100.




Images


Western Blot: KIAA0020 Antibody [NBP1-87296] - Analysis in control (vector only transfected HEK293T lysate) and KIAA0020 over-expression lysate (Co-expressed with a C-terminal myc-DDK tag (3.1 kDa) in mammalian HEK293T cells).




Immunocytochemistry/Immunofluorescence: KIAA0020 Antibody [NBP1-87296] - Immunofluorescent staining of human cell line A-431 shows localization to nucleoli. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: KIAA0020 Antibody [NBP1-87296] - Staining of human colon shows weak to moderate positivity in nucleoli in glandular cells.




Western Blot: KIAA0020 Antibody [NBP1-87296] - Analysis in human cell line HDLM-2.




Immunohistochemistry-Paraffin: KIAA0020 Antibody [NBP1-87296] - Staining of human breast shows nuclear membranous positivity in glandular cells.



Immunohistochemistry-Paraffin: KIAA0020 Antibody [NBP1-87296] - Staining of human kidney shows weak to moderate positivity in nucleoli in cells in tubules.



Immunohistochemistry-Paraffin: KIAA0020 Antibody [NBP1-87296] - Staining of human testis shows weak to moderate positivity in nucleoli in cells in seminiferous ducts.



Immunohistochemistry-Paraffin: KIAA0020 Antibody [NBP1-87296] - Staining of human fallopian tube shows weak to moderate positivity in nucleoli in glandular cells.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-87296

NBP1-87296PEP	PUM3 Recombinant Protein Antigen
NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-87296

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

