

Product Datasheet

Serpin B5/Maspin Antibody - BSA Free NBP1-87778

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-87778

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-87778



NBP1-87778

Serpine B5/Maspin Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Description	Novus Biologicals Rabbit Serpin B5/Maspin Antibody - BSA Free (NBP1-87778) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5268
Gene Symbol	SERPINB5
Species	Human, Rat
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: Mouse (89%).
Immunogen	This antibody was developed against Recombinant Protein corresponding to amino acids: KSLNLSTEFISSTKRPYAKELETVDKDKLEETKGQINNSIKDLTDGHFENILADN SVNDQTKILVVNAAYFVGKWMKKFPESETKECPFRVNKT


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Knockdown Validated
Recommended Dilutions	Western Blot 0.04-0.4 µg/ml, Immunohistochemistry 1:200 - 1:500, Immunohistochemistry-Paraffin 1:200 - 1:500, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Western Blot: Serpin B5/Maspin Antibody [NBP1-87778] - Western blot shows lysates of HeLa human cervical epithelial carcinoma parental cell line and MASPIN/Serpin B5 knockout (KO) HeLa cell line. PVDF membrane was probed with 0.5 µg/ml of Rabbit Anti-Human MASPIN/Serpin B5 Polyclonal Antibody (Catalog # NBP1-87778) followed by HRP-conjugated Anti-Rabbit IgG Secondary Antibody (Catalog #HAF008). Specific band was detected for MASPIN/Serpin B5 at approximately 40 kDa (as indicated) in the parental HeLa cell line, but is not detectable in the knockout HeLa cell line. This experiment was conducted under reducing conditions.




Western Blot: Serpin B5/Maspin Antibody [NBP1-87778] - Analysis in HeLa cells transfected with control siRNA, target specific siRNA probe #1 and #2,. Remaining relative intensity is presented. Loading control: Anti-PPIB.




Immunohistochemistry-Paraffin: Serpin B5/Maspin Antibody [NBP1-87778] - Staining of human liver.




Western Blot: Serpin B5/Maspin Antibody [NBP1-87778] - Analysis in mouse cell line NIH-3T3 and rat cell line NBT-II.




Western Blot: Serpin B5/Maspin Antibody [NBP1-87778] - Analysis in human cell lines PC-3 and Caco-2 using Anti-SERPINB5 antibody. Corresponding SERPINB5 RNA-seq data are presented for the same cell lines. Loading control: Anti-HSP90B1.




Immunohistochemistry-Paraffin: Serpin B5/Maspin Antibody [NBP1-87778] - Staining of human cerebral cortex, liver, skin and testis using Anti-SERPINB5 antibody NBP1-87778 (A) shows similar protein distribution across tissues to independent antibody NBP1-87780 (B).




Immunohistochemistry-Paraffin: Serpin B5/Maspin Antibody [NBP1-87778] - Staining of human skin.



Immunohistochemistry-Paraffin: Serpin B5/Maspin Antibody [NBP1-87778] - Staining of human cerebral cortex.



Immunohistochemistry-Paraffin: Serpin B5/Maspin Antibody [NBP1-87778] - Staining of human testis.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-87778

NBP1-87778PEP	Serpin B5/Maspin Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-87778

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

