

Product Datasheet

Human CD31/PECAM-1 ELISA Kit (Colorimetric) NBP1-91175

Unit Size: 96 Tests

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-91175

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-91175



NBP1-91175**Human CD31/PECAM-1 ELISA Kit (Colorimetric)**


Product Information	
Unit Size	96 Tests
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.

Product Description	
Description	The human sCD31/PECAM-1 ELISA is an enzyme-linked immunosorbent assay for the quantitative detection of human sCD31/PECAM-1
Gene ID	5175
Gene Symbol	PECAM1
Species	Human
Specificity/Sensitivity	The interference of circulating factors of the immune system was evaluated by spiking these proteins at physiologically relevant concentrations into a human CD31/PECAM-1 positive serum. There was no crossreactivity detected.
Standard Curve Range	0.47 - 30.0 ng/mL
Sensitivity	0.06 ng/mL
Inter Assay	CV% < 7.40%
Intra Assay	CV% < 1.70%
Assay Type	Sandwich ELISA
Spiking Recovery	The recovery ranged from 100% to 119% with an overall mean recovery of 108%. Please see protocol for more specific Spike Recovery details.
Suitable Sample Type	Cell culture supernatant, serum, plasma (EDTA, citrate, heparin), urine, other body fluids
Sample Volume	10 uL

Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA

Images

N/A: Human, CD31/PECAM-1 ELISA Kit (Colorimetric) [NBP1-91175] - Standard curve for human Human CD31/PECAM-1 ELISA Kit (Colorimetric)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-91175

NBP2-54655PEP	CD31/PECAM-1 Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]
AF3628	CD31/PECAM-1 Antibody [Unconjugated]
DVE00	VEGF [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-91175

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

