

Product Datasheet

Doublecortin Antibody (3E1) - BSA Free NBP1-92684

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-92684

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-92684




NBP1-92684

Doublecortin Antibody (3E1) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3E1
Preservative	0.035% Sodium Azide
Isotype	IgG2a
Purity	Protein G purified
Buffer	50% PBS, 50% glycerol
Target Molecular Weight	35-45 kDa
Product Description	
Host	Mouse
Gene ID	1641
Gene Symbol	DCX
Species	Human, Mouse, Rat, Porcine, Bovine, Equine
Marker	Immature Neuronal Marker
Immunogen	Recombinant full length human Lis-A isoform of Doublecortin purified from E. coli. [UniProt# O43602]
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:10000, Immunohistochemistry 1:1000, Immunocytochemistry/ Immunofluorescence 1:1000
Application Notes	This Doublecortin (3E1) antibody is useful for Immunocytochemistry/Immunofluorescence, Immunohistochemistry, and Western blot, where bands can be seen at approximately 35 and 45 kDa.

Images


Western Blot: Doublecortin Antibody (3E1) [NBP1-92684] - Analysis of rat whole brain lysates using mouse mAb to Doublecortin (DCX), NBP1-92684, dilution 1:1,000 in green: [1] protein standard (red), [2] adult rat brain, [3] embryonic E20 rat brain. Strong bands at 40, 45kDa correspond to DCX protein, detected only in the developing brain.



Immunocytochemistry/Immunofluorescence: Doublecortin Antibody (3E1) [NBP1-92684] - Analysis of cortical neuron-glia cell culture from E20 rat stained with mouse mAb to Doublecortin (DCX), NBP1-92684, dilution 1:1,000 in red, and costained with chicken pAb to Microtubule Associated Protein2 (MAP2), dilution 1:10,000 in green. The blue is DAPI staining of nuclear DNA. DCX antibody reveals strong cytoplasmic staining in a population of small developing neurons and their processes, while MAP2 antibody stains dendrites and perikarya of mature neurons. Therefore DCX antibody is an excellent marker of developing neuronal cells.



Immunocytochemistry/Immunofluorescence: Doublecortin Antibody (3E1) [NBP1-92684] - Rat brain neural cultures stained with NBP1-92684 (green), a chicken polyclonal antibody to GFAP (NB300-213, red) and DNA (blue). The NBP1-92684 antibody reveals strong cytoplasmic staining in a population of small developing neurons. These cells are often found in small clumps as in these cultures, as in this example. Note that they are not positive for MAP2, which is characteristic of more mature neurons, and that the mature neuron shown is negative for NBP1-92684. This Doublecortin antibody is therefore an excellent marker of developing neuronal cells.



Publications

Obara K, Shirai K, Hamada Y et al. Direct implantation of hair-follicle-associated pluripotent (HAP) stem cells repairs intracerebral hemorrhage and reduces neuroinflammation in mouse model PLoS one 2023-01-13 [PMID: 36638123] (ICC/IF, Mouse)

Details:

Dilution used in IHC-P 1:1000

Kubis-Kubiak A, Wiatrak B, Piwowar A. Et al. Hyper-glycemia and insulinemia induce morphological changes and modulate secretion of S100B, S100A8, amyloid b 1-40 and amyloid b 1-42, in a model of human dopaminergic neurons Biomed Pharmacother 2022-10-18 [PMID: 36257211] (ICC/IF, Rat)

Details:

Citation using the DyLight 488 version of this antibody.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-92684

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)
NBP1-92684G	Doublecortin Antibody (3E1) [DyLight 488]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-92684

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

