

Product Datasheet

Park7/DJ-1 Antibody (4H4) - BSA Free NBP1-92715

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-92715

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-92715




NBP1-92715

Park7/DJ-1 Antibody (4H4) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	4H4
Preservative	0.035% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	50% PBS, 50% glycerol
Target Molecular Weight	21 kDa
Product Description	
Host	Mouse
Gene ID	11315
Gene Symbol	PARK7
Species	Human, Bovine, Mouse (Negative), Rat (Negative)
Reactivity Notes	Not reactive with rodent.
Immunogen	Full length recombinant human Park7(DJ-1) expressed in and purified from E. coli. [UniProt# Q99497]
Product Application Details	
Applications	Western Blot, Simple Western, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Knockout Validated
Recommended Dilutions	Western Blot 1:5000, Simple Western 1:500, Immunohistochemistry 1:5000, Immunocytochemistry/ Immunofluorescence 1:5000, Knockout Validated
Application Notes	In Simple Western only 10 - 15 uL of the recommended dilution is used per data point. See Simple Western Antibody Database for Simple Western validation: Tested in HeLa lysate 0.5 mg/mL, separated by Size, antibody dilution of 1:500, apparent MW was 26 kDa. Separated by Size-Wes, Sally Sue/Peggy Sue.

Images


Western Blot: Park7/DJ-1 Antibody (4H4) [NBP1-92715] - Western blot shows lysates of 293T human embryonic kidney parental cell line and PARK7 knockout (KO) 293T cell line. PVDF membrane was probed with Mouse Anti-Human PARK7 Monoclonal Antibody (Catalog # NBP1-92715) followed by HRP-conjugated Anti-Rabbit IgG Secondary Antibody (Catalog #HAF008). Specific band was detected for PARK7 at approximately 23 kDa (as indicated) in the parental 293T cell line, but is not detectable in the knockout 293T cell line. This experiment was conducted under reducing conditions.



Western Blot: Park7/DJ-1 Antibody (4H4) [NBP1-92715] - Analysis of whole brain and cell lysates using mouse mAb against DJ1/Park7, NBP1-92715, dilution 1:5,000 in green. [1] protein standard, [2] rat brain, [3] mouse brain, [4] NIH-3T3, [5] HeLa, [6] HEK293, and [7] C6 cells. The DJ1 antibody detects protein with apparent molecular weight of 21kDa but only in human cell lines, since it does not recognize the rat or mouse DJ1 protein. The blot was simultaneously probed with chicken pAb to vimentin, dilution 1:5,000 in red, revealing a single band at about 50kDa present in all lanes, though at much lower levels in the tissue lysates.



Immunocytochemistry/Immunofluorescence: Park7(DJ-1) Antibody (4H4) [NBP1-92715] - HeLa cells stained with NBP1-92715 (green), chicken antibody to Vimentin (NB300-223, red) and DNA (blue). NBP1-92715 reveals strong cytoplasmic staining for Park7(DJ-1).



Simple Western: Park7/DJ-1 Antibody (4H4) [NBP1-92715] - Lane view shows a specific band for Park7 (DJ-1) in 0.5 mg/ml of HeLa lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system.



Publications

Repici M, Hassanjani M, Maddison DC et al. The Parkinson's Disease-Linked Protein DJ-1 Associates with Cytoplasmic mRNP Granules During Stress and Neurodegeneration. *Mol Neurobiol* 2018-04-19 [PMID: 29675578]

Wang C, Fang M, Zhang M et al. The positive correlation between DJ-1 and beta-catenin expression shows prognostic value for patients with glioma. *Neuropathology* 2013-05-28 [PMID: 23714193] (ICC/IF, Human)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-92715

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP1-92715F	Park7/DJ-1 Antibody (4H4) [FITC]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-92715

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

