

Product Datasheet

Sheep IgG F(ab)2 Isotype Control [HRP] NBP1-97124

Unit Size: 1 mg

Store lyophilized antibody at 4C in the dark. Aliquot reconstituted liquid and store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-97124

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-97124



NBP1-97124**Sheep IgG F(ab)₂ Isotype Control [HRP]**

| Product Information | |
|------------------------------------|---|
| Unit Size | 1 mg |
| Concentration | LYOPH mg/ml |
| Storage | Store lyophilized antibody at 4C in the dark. Aliquot reconstituted liquid and store at -20C. Avoid freeze-thaw cycles. |
| Preservative | 0.01% Gentamicin Sulfate |
| Reconstitution Instructions | Reconstitute with 1.0 ml deionized water (or equivalent). |
| Isotype | IgG |
| Conjugate | HRP |
| Purity | Multi-step |
| Buffer | Lyophilized from 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free |

| Product Description | |
|--------------------------------|--|
| Description | <p>Store vial at 4C prior to restoration. For extended storage aliquot contents and freeze at -20C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4C as an undiluted liquid. Dilute only prior to immediate use.</p> <p>This product was prepared from normal serum by delipidation, salt fractionation, ion exchange chromatography followed by pepsin digestion and extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Sheep IgG, anti-Sheep IgG F(ab')₂ and anti-Sheep Serum</p> |
| Host | Sheep |
| Species | Sheep |
| Specificity/Sensitivity | No reaction was observed against anti-Sheep IgG F(c) or anti-Pepsin. |

| Product Application Details | |
|------------------------------------|---|
| Applications | Western Blot, ELISA, Immunohistochemistry |
| Recommended Dilutions | Western Blot, ELISA, Immunohistochemistry |
| Application Notes | <p>This product reagents are ideal for ELISA, western blotting, Immunohistochemistry, as well as other antibody detection methods.</p> <p>NOTE: Do NOT add Sodium Azide (it inhibits HRP irreversibly).</p> |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Isotype Controls are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-97124

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

